

BİOSİMA;

enzimler, antijenler,
antikorlar, hızlı testler
ve diğer fonksiyonel
proteinler ve polimer
organik materyallere
yönelik teknoloji ve
ürünlerin ithalatına
odaklanmış bir
biyoteknoloji firmasıdır.

2019 yılında pandeminin
başlaması ile birlikte
medikal ürünlerin ithalat
ve ihracatına odaklanan
firma 2021 yılı itibari ile
Biosima ismi ile
piyasada yerini almıştır.

CATALOG



BİOSİMA



Product Series

Product Contents



Pipette Tips
Automatic Pipette Tips & Universal Pipette Tips



Plate
Deep Well Plate / PCR Plate / Silicone Mat / Sealing Film



Reservoir



Tube/Column



Reagent Bottle



Serological Pipettes



Frozen Storage

Pipette Tips



SBS International Standard

All products conform to SBS international standard, are compatible with most of pipette workstations.



Imported Virgin Material

Imported medical grade polypropylene (PP).



High Temperature Sterilization

Autoclavable at 121°C/15psi.



Low Retention

Smooth tips, no bending, low adsorption, low residue.



Non - Pyrogenic

Non - pyrogen, DNase/RNase Free.



Chemical Resistance

Resistant to chemical reagent during use and storage.



Low Heavy Metal Content

High quality and low heavy metal content.



Traceable

All product package is coded for traceability.

The Applications Of Universal Pipette Tips

> DNA Extraction

DNA extraction (e.g. PCR).

> RNA Applications

RNA applications (e.g. gene analysis).

> Protein Applications

Protein applications (e.g. antibody studies).

> Cell Culture Applications

Cell culture applications (e.g. media).

> Aerosol Bonding Liquid

All applications using aerosol-bonded liquids.



Pipette Tip Package



Individually Rack

Easy to use

96 pcs/rack



Stacked Box Packaging

Save space and reduce plastic waste by 60%

96 pcs/layer



Bagged Packing

Save space and reduce plastic waste by 85%

1000 pcs/bag

Variety of sizes, types and packaging are available to meet needed.

■ 10 μ L / 10 μ L Extra Long Tips

■ 200 μ L / 300 μ L / 300 μ L Extra Long Tips

■ 1000 μ L / 1250 μ L Tips




Compatibility List For Agilent And Amtk Tips



AMTK



Package: 384 pcs/rack , 96 pcs/rack
Volume: 30uL , 70uL , 250uL
Material: Medical grade polypropylene
Temperature range: -80°C-121°C



Agilent Pipette Tips	Biosima Model	Volume(μL)	Sterile	Low-Retention	Filtering	Packing
	VT-384-30-RS	30	✓			384 pcs/rack, 10 racks /box, 5 boxes/carton
	VT-384-30-RSL	30	✓	✓		384 pcs/rack, 10 racks /box, 5 boxes/carton
	VT-384-70-RS	70	✓			384 pcs/rack, 10 racks /box, 5 boxes/carton
	VT-384-70-RSL	70	✓	✓		384 pcs/rack, 10 racks /box, 5 boxes/carton
	VT-250-R	250				96 pcs/rack, 10 racks /box, 5 boxes/carton
	VT-250-RL	250		✓		96 pcs/rack, 10 racks /box, 5 boxes/carton
	VT-250-RS	250	✓			96 pcs/rack, 10 racks /box, 5 boxes/carton
	VT-250-RSL	250	✓	✓		96 pcs/rack, 10 racks /box, 5 boxes/carton
	VTF-250-RS	250	✓		✓	96 pcs/rack, 10 racks /box, 5 boxes/carton
	VTF-250-RSL	250	✓	✓	✓	96 pcs/rack, 10 racks /box, 5 boxes/carton

Compatible with Agilent Bravo and Vprep / AMTK , LH-1510/LH-1510+ automatic liquid handling workstation.

Compatibility List For Beckman And Nayo Tips



Package: 96 pcs/rack
Volume: 50uL , 250uL
Material: Medical grade polypropylene
Temperature range: -80°C-121°C

Beckman Pipette Tips	Biosima Model	Volume(μL)	Sterile	Low-Retention	Filtering	Packing
	BK-50-R	50				96 pcs/rack, 10 racks /box, 5 boxes/carton
	BK-50-RL	50		✓		96 pcs/rack, 10 racks /box, 5 boxes/carton
	BK-50-RS	50	✓			96 pcs/rack, 10 racks /box, 5 boxes/carton
	BK-50-RSL	50	✓	✓		96 pcs/rack, 10 racks /box, 5 boxes/carton
	BKF-50-RS	50	✓		✓	96 pcs/rack, 10 racks /box, 5 boxes/carton
	BKF-50-RSL	50	✓	✓	✓	96 pcs/rack, 10 racks /box, 5 boxes/carton
	BK-250-R	250				96 pcs/rack, 10 racks /box, 5 boxes/carton
	BK-250-RL	250		✓		96 pcs/rack, 10 racks /box, 5 boxes/carton
	BK-250-RS	250	✓			96 pcs/rack, 10 racks /box, 5 boxes/carton
	BK-250-RSL	250	✓	✓		96 pcs/rack, 10 racks /box, 5 boxes/carton
	BKF-250-RS	250	✓		✓	96 pcs/rack, 10 racks /box, 5 boxes/carton
	BKF-250-RSL	250	✓	✓	✓	96 pcs/rack, 10 racks /box, 5 boxes/carton

Compatible with Biomek i5 , Biomek i7 NAYO , NAYO 96 automatic liquid handling workstation.

Compatible List For Tecan Transparent Tips





TECAN.




ALLSHENG



Package: 96 pcs/rack
Volume: 20uL , 50uL , 200uL , 1000uL
Material: Medical grade polypropylene
Temperature range: -80°C-121°C

Tecan Tips (Transparent)	Biosima Model	Volume(μL)	Sterile	Low-Retention	Filtering	Packing
	TT-20-R	20				96 pcs/rack, 24 racks/carton
	TT-20-RL	20		✓		96 pcs/rack, 24 racks/carton
	TT-20-RS	20	✓			96 pcs/rack, 24 racks/carton
	TT-20-RSL	20	✓	✓		96 pcs/rack, 24 racks/carton
	TTF-20-RS	20	✓		✓	96 pcs/rack, 24 racks/carton
	TTF-20-RSL	20	✓	✓	✓	96 pcs/rack, 24 racks/carton
	TT-50-R	50				96 pcs/rack, 24 racks/carton
	TT-50-RL	50		✓		96 pcs/rack, 24 racks/carton
	TT-50-RS	50	✓			96 pcs/rack, 24 racks/carton
	TT-50-RSL	50	✓	✓		96 pcs/rack, 24 racks/carton
	TTF-50-RS	50	✓		✓	96 pcs/rack, 24 racks/carton
	TTF-50-RSL	50	✓	✓	✓	96 pcs/rack, 24 racks/carton


Compatible with Tecan LiHa / ALLSHENG, Au-mate 96 automated pipetting workstation.

Tecan Tips (Transparent)	Biosima Model	Volume(μL)	Sterile	Low-Retention	Filtering	Packing
	TT-200-R	200				96 pcs/rack, 24 racks/carton
	TT-200-RL	200		✓		96 pcs/rack, 24 racks/carton
	TT-200-RS	200	✓			96 pcs/rack, 24 racks/carton
	TT-200-RSL	200	✓	✓		96 pcs/rack, 24 racks/carton
	TTF-200-RS	200	✓		✓	96 pcs/rack, 24 racks/carton
	TTF-200-RSL	200	✓	✓	✓	96 pcs/rack, 24 racks/carton

Compatible with Tecan LiHa / ALLSHENG, Au-mate 96 automated pipetting workstation.

- 1 Tecan Tips (20uL)
Type: sterilized,low retention,Filtering
- 2 Tecan Tips (50uL)
Type: sterilized,low retention,Filtering
- 3 Tecan Tips (200uL)
Type: sterilized,low retention,Filtering
- 4 Tecan Tips (1000uL)
Type: sterilized,low retention,Filtering



Tecan Tips (Transparent)	Biosima Model	Volume(μL)	Sterile	Low-Retention	Filtering	Packing
	TT-1000-R	1 000				96 pcs/rack, 16 racks/carton
	TT-1000-RL	1 000		✓		96 pcs/rack, 16 racks/carton
	TT-1000-RS	1 000	✓			96 pcs/rack, 16 racks/carton
	TT-1000-RSL	1 000	✓	✓		96 pcs/rack, 16 racks/carton
	TTF-1000-RS	1 000	✓		✓	96 pcs/rack, 16 racks/carton
	TTF-1000-RSL	1 000	✓	✓	✓	96 pcs/rack, 16 racks/carton

Compatible with Tecan LiHa / ALLSHENG, Au-mate 96 automated pipetting workstation.

Compatible List For Tecan Transparent Tips (SBS Rack Box)





TECAN.




ALLSHENG



Package: 96 pcs/rack
Volume: 20uL , 50uL , 200uL , 1000uL
Material: Medical grade polypropylene
Temperature range: -80°C-121°C

Tecan Tips (Transparent)	Biosima Model	Volume(μL)	Sterile	Low-Retention	Filtering	Packing
	TT-20-H	20				96 pcs/rack, 24 racks/carton
	TT-20-HL	20		✓		96 pcs/rack, 24 racks/carton
	TT-20-HS	20	✓			96 pcs/rack, 24 racks/carton
	TT-20-HSL	20	✓	✓		96 pcs/rack, 24 racks/carton
	TTF-20-HS	20	✓		✓	96 pcs/rack, 24 racks/carton
	TTF-20-HSL	20	✓	✓	✓	96 pcs/rack, 24 racks/carton
	TT-50-H	50				96 pcs/rack, 24 racks/carton
	TT-50-HL	50		✓		96 pcs/rack, 24 racks/carton
	TT-50-HS	50	✓			96 pcs/rack, 24 racks/carton
	TT-50-HSL	50	✓	✓		96 pcs/rack, 24 racks/carton
	TTF-50-HS	50	✓		✓	96 pcs/rack, 24 racks/carton
	TTF-50-HSL	50	✓	✓	✓	96 pcs/rack, 24 racks/carton

Compatible with Tecan LiHa / ALLSHENG, Auto-Pure 4800 automated pipetting workstation.

Tecan Tips (Transparent)	Biosima Model	Volume(μL)	Sterile	Low-Retention	Filtering	Packing
	TT-200-H	200				96 pcs/rack, 24 racks/carton
	TT-200-HL	200		✓		96 pcs/rack, 24 racks/carton
	TT-200-HS	200	✓			96 pcs/rack, 24 racks/carton
	TT-200-HSL	200	✓	✓		96 pcs/rack, 24 racks/carton
	TTF-200-HS	200	✓		✓	96 pcs/rack, 24 racks/carton
	TTF-200-HSL	200	✓	✓	✓	96 pcs/rack, 24 racks/carton

Compatible with Tecan LiHa / ALLSHENG, Auto-Pure 4800 automated pipetting workstation.


- 1


Tecan Tips (20uL)
Type: sterilized,low retention,Filtering
- 2

Tecan Tips (50uL)
Type: sterilized,low retention,Filtering
- 3

Tecan Tips (200uL)
Type: sterilized,low retention,Filtering
- 4

Tecan Tips (1000uL)
Type: sterilized,low retention,Filtering



Tecan Tips (Transparent)	Biosima Model	Volume(μL)	Sterile	Low-Retention	Filtering	Packing
	TT-1000-H	1 000				96 pcs/rack, 16 racks/carton
	TT-1000-HL	1 000		✓		96 pcs/rack, 16 racks/carton
	TT-1000-HS	1 000	✓			96 pcs/rack, 16 racks/carton
	TT-1000-HSL	1 000	✓	✓		96 pcs/rack, 16 racks/carton
	TTF-1000-HS	1 000	✓		✓	96 pcs/rack, 16 racks/carton
	TTF-1000-HSL	1 000	✓	✓	✓	96 pcs/rack, 16 racks/carton

Compatible with Tecan LiHa / ALLSHENG, Auto-Pure 4800 automated pipetting workstation.

Compatible List For Tecan Conductive Tips





TECAN.




ALLSHENG



Package: 96 pcs/rack
Volume: 20uL , 50uL , 200uL , 1000uL
Material: Medical grade polypropylene
Temperature range: -80°C-121°C


Tecan Tips (conductive)	Biosima Model	Volume(μL)	Sterile	Low-Retention	Filtering	Packing
	TT-20C-R	20				96 pcs/rack, 24 racks/carton
	TT-20C-RL	20		✓		96 pcs/rack, 24 racks/carton
	TT-20C-RS	20	✓			96 pcs/rack, 24 racks/carton
	TT-20C-RSL	20	✓	✓		96 pcs/rack, 24 racks/carton
	TTF-20C-RS	20	✓		✓	96 pcs/rack, 24 racks/carton
	TTF-20C-RSL	20	✓	✓	✓	96 pcs/rack, 24 racks/carton
	TT-50C-R	50				96 pcs/rack,24 racks/carton
	TT-50C-RL	50		✓		96 pcs/rack,24 racks/carton
	TT-50C-RS	50	✓			96 pcs/rack,24 racks/carton
	TT-50C-RSL	50	✓	✓		96 pcs/rack, 24 racks/carton
	TTF-50C-RS	50	✓		✓	96 pcs/rack, 24 racks/carton
	TTF-50C-RSL	50	✓	✓	✓	96 pcs/rack, 24 racks/carton

Compatible with Tecan LiHa, For Genesis Freedom, Freedom Evo / ALLSHENG, Auto-Pure 4800 automated pipetting workstation.

Tecan Tips (conductive)	Biosima Model	Volume(μL)	Sterile	Low-Retention	Filtering	Packing
	TT-200C-R	200				96 pcs/rack, 24 racks/carton
	TT-200C-RL	200		✓		96 pcs/rack, 24 racks/carton
	TT-200C-RS	200	✓			96 pcs/rack, 24 racks/carton
	TT-200C-RSL	200	✓	✓		96 pcs/rack, 24 racks/carton
	TTF-200C-RS	200	✓		✓	96 pcs/rack, 24 racks/carton
	TTF-200C-RSL	200	✓	✓	✓	96 pcs/rack, 24 racks/carton

Compatible with Tecan LiHa, For Genesis Freedom, Freedom Evo / ALLSHENG, Auto-Pure 4800 automated pipetting workstation.



Tecan Tips (conductive)	Biosima Model	Volume(μL)	Sterile	Low-Retention	Filtering	Packing
	TT-1000C-R	1 000				96 pcs/rack, 16 racks/carton
	TT-1000C-RL	1 000		✓		96 pcs/rack, 16 racks/carton
	TT-1000C-RS	1 000	✓			96 pcs/rack, 16 racks/carton
	TT-1 000C-RSL	1 000	✓	✓		96 pcs/rack, 16 racks/carton
	TTF-1000C-RS	1 000	✓		✓	96 pcs/rack, 16 racks/carton
	TTF-1000C-RSL	1 000	✓	✓	✓	96 pcs/rack, 16 racks/carton

Compatible with Tecan LiHa, For Genesis Freedom, Freedom Evo / ALLSHENG, Auto-Pure 4800 automated pipetting workstation.

Compatible List For Tecan
Conductive Tips (SBS Rack Box)





TECAN.





ALLSHENG



Package: 96 pcs/rack
Volume: 20uL , 50uL , 200uL ,
1000uL , 5000uL
Material: Medical grade
polypropylene
Temperature range: -80°C-121°C


Tecan Tips (conductive)	Biosima Model	Volume(μL)	Sterile	Low-Retention	Filtering	Packing
	TT-20C-H	20				96 pcs/rack, 24 racks/carton
	TT-20C-HL	20		✓		96 pcs/rack, 24 racks/carton
	TT-20C-HS	20	✓			96 pcs/rack, 24 racks/carton
	TT-20C-HSL	20	✓	✓		96 pcs/rack, 24 racks/carton
	TTF-20C-HS	20	✓		✓	96 pcs/rack, 24 racks/carton
	TTF-20C-HSL	20	✓	✓	✓	96 pcs/rack, 24 racks/carton
	TT-50C-H	50				96 pcs/rack, 24 racks/carton
	TT-50C-HL	50		✓		96 pcs/rack, 24 racks/carton
	TT-50C-HS	50	✓			96 pcs/rack, 24 racks/carton
	TT-50C-HSL	50	✓	✓		96 pcs/rack, 24 racks/carton
	TTF-50C-HS	50	✓		✓	96 pcs/rack, 24 racks/carton
	TTF-50C-HSL	50	✓	✓	✓	96 pcs/rack, 24 racks/carton

Compatible with Tecan LiHa, For Genesis Freedom, Freedom Evo / ALLSHENG, Auto-Pure 4800 automated pipetting workstation.

Tecan Tips (conductive)	Biosima Model	Volume(μL)	Sterile	Low-Retention	Filtering	Packing
	TT-200C-H	200				96 pcs/rack, 24 racks/carton
	TT-200C-HL	200		✓		96 pcs/rack, 24 racks/carton
	TT-200C-HS	200	✓			96 pcs/rack, 24 racks/carton
	TT-200C-HSL	200	✓	✓		96 pcs/rack, 24 racks/carton
	TTF-200C-HS	200	✓		✓	96 pcs/rack, 24 racks/carton
	TTF-200C-HSL	200	✓	✓	✓	96 pcs/rack, 24 racks/carton
	TT-1000C-H	1000				96 pcs/rack, 16 racks/carton
	TT-1000C-HL	1000		✓		96 pcs/rack, 16 racks/carton
	TT-1000C-HS	1000	✓			96 pcs/rack, 16 racks/carton
	TT-1000C-HSL	1000	✓	✓		96 pcs/rack, 16 racks/carton
	TTF-1000C-HS	1000	✓		✓	96 pcs/rack, 16 racks/carton
	TTF-1000C-HSL	1000	✓	✓	✓	96 pcs/rack, 16 racks/carton


- 1 Tecan Tips (20uL)
Type: sterilized,low retention,Filtering
- 2 Tecan Tips (50uL)
Type: sterilized,low retention,Filtering
- 3 Tecan Tips (200uL)
Type: sterilized,low retention,Filtering
- 4 Tecan Tips (1000uL)
Type: sterilized,low retention,Filtering
- 5 Tecan Tips (5000uL)
Type: sterilized,low retention,Filtering



Tecan Tips (conductive)	Biosima Model	Volume(μL)	Sterile	Retention	Filtering	Packing
	TT-5000C-H	5000				96 pcs/rack, 16 racks/carton
	TT-5000C-HL	5000		✓		96 pcs/rack, 16 racks/carton
	TT-5000C-HS	5000	✓			96 pcs/rack, 16 racks/carton
	TT-5000C-HSL	5000	✓	✓		96 pcs/rack, 16 racks/carton
	TTF-5000C-HS	5000	✓		✓	96 pcs/rack, 16 racks/carton
	TTF-5000C-HSL	5000	✓	✓	✓	96 pcs/rack, 16 racks/carton

Compatible with Tecan LiHa, For Genesis Freedom, Freedom Evo / ALLSHENG, Auto-Pure 4800 automated pipetting workstation.



Compatibility List For Hamilton
Conductive Tips



HAMILTON



Package: 96 pcs/rack
Volume: 50uL , 300uL , 1000uL
Material: Medical grade polypropylene
Temperature range: -80°C-121°C


Hamilton Tips (conductive)	Biosima Model	Volume(μL)	Sterile	Low-Retention	Filtering	Packing
	HT-50C-RS	50	✓			96 pcs/rack, 36 racks/carton
	HT-50C-RSL	50	✓	✓		96 pcs/rack, 36 racks/carton
	HTF-50C-RS	50	✓		✓	96 pcs/rack, 36 racks/carton
	HTF-50C-RSL	50	✓	✓	✓	96 pcs/rack, 36 racks/carton
	HT-300C-RS	300	✓			96 pcs/rack, 36 racks/carton
	HT-300C-RSL	300	✓	✓		96 pcs/rack, 36 racks/carton
	HTF-300C-RS	300	✓		✓	96 pcs/rack, 36 racks/carton
	HTF-300C-RSL	300	✓	✓	✓	96 pcs/rack, 36 racks/carton
	HT-1000C-RS	1000	✓			96 pcs/rack, 24 racks/carton
	HT-1000C-RSL	1000	✓	✓		96 pcs/rack, 24 racks/carton
	HTF-1000C-RS	1000	✓		✓	96 pcs/rack, 24 racks/carton
	HTF-1000C-RSL	1000	✓	✓	✓	96 pcs/rack, 24 racks/carton

Compatible with the liquid handling workstation of Hamilton Star, Starlet and N imbus.




Compatibility List For
Opentrons Tips




opentrons






Package: 96 pcs/rack
Volume: 20uL , 200uL , 1000uL
Material: Medical grade polypropylene
Temperature range: -80°C-121°C





Opentrons Pipette Tips	Biosima Model	Volume(μL)	Sterile	Low-Retention	Filtering	Packing
	OPT-20-RS	20	✓			96 pcs/rack, 10 racks /box, 5 boxes /carton
	OPT-20-RSL	20	✓	✓		96 pcs/rack, 10 racks /box, 5 boxes /carton
	OPTF-20-RS	20	✓		✓	96 pcs/rack, 10 racks /box, 5 boxes /carton
	OPTF-20-RSL	20	✓	✓	✓	96 pcs/rack, 10 racks /box, 5 boxes /carton
	OPT-20-ES	20	✓			96 pcs /layer, 10 layers /box (reloading package)
	OPT-20-ESL	20	✓	✓		96 pcs /layer, 10 layers /box (reloading package)
	OPT-200-RS	200	✓			96 pcs/rack, 10 racks /box, 5 boxes /carton
	OPT-200-RSL	200	✓	✓		96 pcs/rack, 10 racks /box, 5 boxes /carton
	OPTF-200-RS	200	✓		✓	96 pcs/rack, 10 racks /box, 5 boxes /carton
	OPTF-200-RSL	200	✓	✓	✓	96 pcs/rack, 10 racks /box, 5 boxes /carton
	OPT-200-ES	200	✓			96 pcs /layer, 10 layers /box (reloading package)
	OPT-200-ESL	200	✓	✓		96 pcs /layer, 10 layers /box (reloading package)
	OPT-1000-RS	1000	✓			96 pcs/rack, 10 racks /box, 5 boxes /carton
	OPT-1000-RSL	1000	✓	✓		96 pcs/rack, 10 racks /box, 5 boxes /carton
	OPTF-1000-RS	1000	✓		✓	96 pcs/rack, 10 racks /box, 5 boxes /carton
	OPTF-1000-RSL	1000	✓	✓	✓	96 pcs/rack, 10 racks /box, 5 boxes /carton
	OPT-1000-ES	1000	✓			96 pcs /layer, 10 layers /box (reloading package)
	OPT-1000-ESL	1000	✓	✓		96 pcs /layer, 10 layers /box (reloading package)



Compatible with the liquid handling workstation of Opentrons, OT-2 Robot.



10μL Tips	Biosima Model	Volume(μL)	Sterile	Low-Retention	Filtering	Packing
 Tip Marks: 3μL, 10μL	T-10-B	0.5-10				1000 pcs/pack, 20 packs/carton
	T-10-BL	0.5-10		✓		1000 pcs/pack, 20 packs/carton
	T-10-R	0.5-10				96 pcs/rack, 10 racks /box, 5 boxes/carton
	T-10-RL	0.5-10		✓		96 pcs/rack, 10 racks /box, 5 boxes/carton
	T-10-RS	0.5-10	✓			96 pcs/rack, 10 racks /box, 5 boxes/carton
	T-10-RSL	0.5-10	✓	✓		96 pcs/rack, 10 racks /box, 5 boxes/carton
	T-10-E	0.5-10				96 pcs/rack, 10 racks /box, 10 boxes/carton
	T-10-EL	0.5-10		✓		96 pcs/rack, 10 racks /box, 10 boxes/carton
10μL Extra Long Tips	Biosima Model	Volume(μL)	Sterile	Low-Retention	Filtering	Packing
 Tip Marks: 2μL, 10μL	T-11-B	0.5-10				1000 pcs/pack, 20 packs/carton
	T-11-BL	0.5-10		✓		1000 pcs/pack, 20 packs/carton
	T-11-R	0.5-10				96 pcs/rack, 10 racks /box, 5 boxes/carton
	T-11-RL	0.5-10		✓		96 pcs/rack, 10 racks /box, 5 boxes/carton
	T-11-RS	0.5-10	✓			96 pcs/rack, 10 racks /box, 5 boxes/carton
	T-11-RSL	0.5-10	✓	✓		96 pcs/rack, 10 racks /box, 5 boxes/carton
10μL Extra Long Tips	Biosima Model	Volume(μL)	Sterile	Low-Retention	Filtering	Packing
 Tip Marks: 2.5μL, 10μL	T-12-B	0.5-10				1000 pcs/pack, 20 packs/carton
	T-12-BL	0.5-10		✓		1000 pcs/pack, 20 packs/carton
	T-12-R	0.5-10				96 pcs/rack, 10 racks /box, 5 boxes/carton
	T-12-RL	0.5-10		✓		96 pcs/rack, 10 racks /box, 5 boxes/carton
	T-12-RS	0.5-10	✓			96 pcs/rack, 10 racks /box, 5 boxes/carton
	T-12-RSL	0.5-10	✓	✓		96 pcs/rack, 10 racks /box, 5 boxes/carton
200μL Tips (Yellow)	Biosima Model	Volume(μL)	Sterile	Low-Retention	Filtering	Packing
	T-200Y-B	10-200				1000 pcs/pack, 20 packs/carton
	T-200Y-BL	10-200		✓		1000 pcs/pack, 20 packs/carton
	T-200Y-R	10-200				96 pcs/rack, 10 racks /box, 5 boxes/carton
	T-200Y-RL	10-200		✓		96 pcs/rack, 10 racks /box, 5 boxes/carton
	T-200Y-RS	10-200	✓			96 pcs/rack, 10 racks /box, 5 boxes/carton
	T-200Y-RSL	10-200	✓	✓		96 pcs/rack, 10 racks /box, 5 boxes/carton
	T-200Y-E	10-200				96 pcs/rack, 10 racks /box, 10 boxes/carton
	T-200Y-EL	10-200		✓		96 pcs/rack, 10 racks /box, 10 boxes/carton



Compatible with a variety of pipettors, suitable for pipetting micro reagent from deep vessels.

300µL Tips	Biosima Model	Volume(µL)	Sterile	Low-Retention	Filtering	Packing
 Tip Marks: 10µL, 50µL, 100µL	T-300-B	10-300				1000 pcs/pack, 20 packs/carton
	T-300-BL	10-300		✓		1000 pcs/pack, 20 packs/carton
	T-300-R	10-300				96 pcs/rack,10 racks /box, 5 boxes/carton
	T-300-RL	10-300		✓		96 pcs/rack,10 racks /box, 5 boxes/carton
	T-300-RS	10-300	✓			96 pcs/rack,10 racks /box, 5 boxes/carton
	T-300-RSL	10-300	✓	✓		96 pcs/rack,10 racks /box, 5 boxes/carton
	T-300-E	10-300				96 pcs/rack,10 racks /box, 10 boxes/carton
	T-300-EL	10-300		✓		96 pcs/rack,10 racks /box, 10 boxes/carton
300µL Extra Long Tips	Biosima Model	Volume(µL)	Sterile	Low-Retention	Filtering	Packing
 Tip Marks: 20µL, 50µL, 100µL, 200µL	T-301-B	10-300				1000 pcs/pack, 5 packs/carton
	T-301-BL	10-300		✓		1000 pcs/pack, 5 packs/carton
	T-301-R	10-300				96 pcs/rack,10 racks /box, 5 boxes/carton
	T-301-RL	10-300		✓		96 pcs/rack,10 racks /box, 5 boxes/carton
	T-301-RS	10-300	✓			96 pcs/rack,10 racks /box, 5 boxes/carton
	T-301-RSL	10-300	✓	✓		96 pcs/rack,10 racks /box, 5 boxes/carton
	T-301-E	10-300				96 pcs/rack, 5 racks /box, 10 boxes/carton
	T-301-EL	10-300		✓		96 pcs/rack, 5 racks /box, 10 boxes/carton
300µL Extra Long Tips	Biosima Model	Volume(µL)	Sterile	Low-Retention	Filtering	Packing
 Tip Marks: 30µL, 100µL, 200µL, 300µL	T-302-B	10-300				1000 pcs/pack, 5 packs/carton
	T-302-BL	10-300		✓		1000 pcs/pack, 5 packs/carton
	T-302-R	10-300				96 pcs/rack,10 racks /box, 5 boxes/carton
	T-302-RL	10-300		✓		96 pcs/rack,10 racks /box, 5 boxes/carton
	T-302-RS	10-300	✓			96 pcs/rack,10 racks /box, 5 boxes/carton
	T-302-RSL	10-300	✓	✓		96 pcs/rack,10 racks /box, 5 boxes/carton
1000µL Tips (Blue)	Biosima Model	Volume(µL)	Sterile	Low-Retention	Filtering	Packing
 Tip Marks: 200µL, 500µL	T-1000B-B	100-1000				1000 pcs/pack, 5 packs/carton
	T-1000B-BL	100-1000		✓		1000 pcs/pack, 5 packs/carton
	T-1000B-R	100-1000				96 pcs/rack,10 racks /box, 5 boxes/carton
	T-1000B-RL	100-1000		✓		96 pcs/rack,10 racks /box, 5 boxes/carton
	T-1000B-RS	100-1000	✓			96 pcs/rack,10 racks /box, 5 boxes/carton
	T-1000B-RSL	100-1000	✓	✓		96 pcs/rack,10 racks /box, 5 boxes/carton
	T-1000B-E	100-1000				96 pcs/rack, 5 racks /box, 10 boxes/carton
	T-1000B-EL	100-1000		✓		96 pcs/rack, 5 racks /box, 10 boxes/carton

1250µL Tips	Biosima Model	Volume(µL)	Sterile	Low-Retention	Filtering	Packing
 Tip Marks: 200µL, 500µL, 1000µL	T-1250-B	100-1250				1000 pcs/pack, 5 packs/carton
	T-1250-BL	100-1250		✓		1000 pcs/pack, 5 packs/carton
	T-1250-R	100-1250				96 pcs/rack,10 racks /box, 5 boxes/carton
	T-1250-RL	100-1250		✓		96 pcs/rack,10 racks /box, 5 boxes/carton
	T-1250-RS	100-1250	✓			96 pcs/rack,10 racks /box, 5 boxes/carton
	T-1250-RSL	100-1250	✓	✓		96 pcs/rack,10 racks /box, 5 boxes/carton
	T-1250-E	100-1250				96 pcs/rack, 5 racks /box, 10 boxes/carton
	T-1250-EL	100-1250		✓		96 pcs/rack, 5 racks /box, 10 boxes/carton
1500µL Tips	Biosima Model	Volume(µL)	Sterile	Low-Retention	Filtering	Packing
 Tip Marks: 500µL, 1000µL, 1500µL	T-1500-B	100-1500				1000 pcs/pack, 5 packs/carton
	T-1500-BL	100-1500		✓		1000 pcs/pack, 5 packs/carton
	T-1500-R	100-1500				96 pcs/rack, 5 racks /box, 5 boxes/carton
	T-1500-RL	100-1500		✓		96 pcs/rack, 5 racks /box, 5 boxes/carton
	T-1500-RS	100-1500	✓			96 pcs/rack, 5 racks /box, 5 boxes/carton
	T-1500-RSL	100-1500	✓	✓		96 pcs/rack, 5 racks /box, 5 boxes/carton

Compatible with a variety of pipettors, suitable for pipetting micro reagent from deep vessels.

Gel Sampling Pipette Tips







Gel Sampling Pipette Tips	Biosima Model	Volume(µL)	Sterile	Low-Retention	Filtering	Packing
	TGL-200-R	200				96 pcs/rack,10 racks /box, 5 boxes/carton
	TGLF-160-RS	160	✓		✓	96 pcs/rack,10 racks /box, 5 boxes/carton

Compatible with a variety of pipettors, suitable for pipetting micro reagent from deep vessels.

PIPETTE TIPS

Universal Pipette Tips
(Filter Element)




10µL Tips	Biosima Model	Volume(µL)	Sterile	Low-Retention	Filtering	Packing
 Tip Marks: 3µL, 10µL	TF-10-B	0.5-10			✓	1000 pcs/pack, 10 packs/carton
	TF-10-BL	0.5-10		✓	✓	1000 pcs/pack, 10 packs/carton
	TF-10-RS	0.5-10	✓		✓	96 pcs/rack, 10 racks /box, 5 boxes/carton
	TF-10-RSL	0.5-10	✓	✓	✓	96 pcs/rack, 10 racks /box, 5 boxes/carton
10µL Extra Long Tips	Biosima Model	Volume(µL)	Sterile	Low-Retention	Filtering	Packing
 Tip Marks: 2µL, 10µL	TF-11-B	0.5-10			✓	1000 pcs/pack, 10 packs/carton
	TF-11-BL	0.5-10		✓	✓	1000 pcs/pack, 10 packs/carton
	TF-11-RS	0.5-10	✓		✓	96 pcs/rack, 10 racks /box, 5 boxes/carton
	TF-11-RSL	0.5-10	✓	✓	✓	96 pcs/rack, 10 racks /box, 5 boxes/carton
10µL Extra Long Tips	Biosima Model	Volume(µL)	Sterile	Low-Retention	Filtering	Packing
 Tip Marks: 2.5µL, 10µL	TF-12-B	0.5-10			✓	1000 pcs/pack, 10 packs/carton
	TF-12-BL	0.5-10		✓	✓	1000 pcs/pack, 10 packs/carton
	TF-12-RS	0.5-10	✓		✓	96 pcs/rack, 10 racks /box, 5 boxes/carton
	TF-12-RSL	0.5-10	✓	✓	✓	96 pcs/rack, 10 racks /box, 5 boxes/carton
20µL Tips (Yellow)	Biosima Model	Volume(µL)	Sterile	Low-Retention	Filtering	Packing
	TF-20Y-B	1.0-20			✓	1000 pcs/pack, 10 packs/carton
	TF-20Y-BL	1.0-20		✓	✓	1000 pcs/pack, 10 packs/carton
	TF-20Y-RS	1.0-20	✓		✓	96 pcs/rack, 10 racks /box, 5 boxes/carton
	TF-20Y-RSL	1.0-20	✓	✓	✓	96 pcs/rack, 10 racks /box, 5 boxes/carton
50µL Tips (Yellow)	Biosima Model	Volume(µL)	Sterile	Low-Retention	Filtering	Packing
	TF-50Y-B	5-50			✓	1000 pcs/pack, 10 packs/carton
	TF-50Y-BL	5-50		✓	✓	1000 pcs/pack, 10 packs/carton
	TF-50Y-RS	5-50	✓		✓	96 pcs/rack, 10 racks /box, 5 boxes/carton
	TF-50Y-RSL	5-50	✓	✓	✓	96 pcs/rack, 10 racks /box, 5 boxes/carton
100µL Tips	Biosima Model	Volume(µL)	Sterile	Low-Retention	Filtering	Packing
	TF-100-B	10-100			✓	1000 pcs/pack, 10 packs/carton
	TF-100-BL	10-100		✓	✓	1000 pcs/pack, 10 packs/carton
	TF-100-RS	10-100	✓		✓	96 pcs/rack, 10 racks /box, 5 boxes/carton
	TF-100-RSL	10-100	✓	✓	✓	96 pcs/rack, 10 racks /box, 5 boxes/carton

Compatible with a variety of pipettors, suitable for pipetting micro reagent from deep vessels.




PIPETTE TIPS

300µL Tips	Biosima Model	Volume(µL)	Sterile	Low-Retention	Filtering	Packing
 Tip Marks: 10µL, 50µL, 100µL	T-300-B	10-300				1000 pcs/pack, 20 packs/carton
	T-300-BL	10-300		✓		1000 pcs/pack, 20 packs/carton
	T-300-R	10-300				96 pcs/rack,10 racks /box, 5 boxes/carton
	T-300-RL	10-300		✓		96 pcs/rack,10 racks /box, 5 boxes/carton
	T-300-RS	10-300	✓			96 pcs/rack,10 racks /box, 5 boxes/carton
	T-300-RSL	10-300	✓	✓		96 pcs/rack,10 racks /box, 5 boxes/carton
	T-300-E	10-300				96 pcs/rack,10 racks/box, 10 boxes/carton
	T-300-EL	10-300		✓		96 pcs/rack,10 racks/box, 10 boxes/carton
300µL Extra Long Tips	Biosima Model	Volume(µL)	Sterile	Low-Retention	Filtering	Packing
 Tip Marks: 20µL, 50µL, 100µL, 200µL	T-301-B	10-300				1000 pcs/pack, 5 packs/carton
	T-301-BL	10-300		✓		1000 pcs/pack, 5 packs/carton
	T-301-R	10-300				96 pcs/rack,10 racks /box, 5 boxes/carton
	T-301-RL	10-300		✓		96 pcs/rack,10 racks /box, 5 boxes/carton
	T-301-RS	10-300	✓			96 pcs/rack,10 racks /box, 5 boxes/carton
	T-301-RSL	10-300	✓	✓		96 pcs/rack,10 racks /box, 5 boxes/carton
	T-301-E	10-300				96 pcs/rack, 5 racks/box, 10 boxes/carton
	T-301-EL	10-300		✓		96 pcs/rack, 5 racks/box, 10 boxes/carton
300µL Extra Long Tips	Biosima Model	Volume(µL)	Sterile	Low-Retention	Filtering	Packing
 Tip Marks: 30µL, 100µL, 200µL, 300µL	T-302-B	10-300				1000 pcs/pack, 5 packs/carton
	T-302-BL	10-300		✓		1000 pcs/pack, 5 packs/carton
	T-302-R	10-300				96 pcs/rack,10 racks /box, 5 boxes/carton
	T-302-RL	10-300		✓		96 pcs/rack,10 racks /box, 5 boxes/carton
	T-302-RS	10-300	✓			96 pcs/rack,10 racks /box, 5 boxes/carton
	T-302-RSL	10-300	✓	✓		96 pcs/rack,10 racks /box, 5 boxes/carton
	T-302-E	10-300				96 pcs/rack, 5 racks/box, 10 boxes/carton
1000µL Tips (Blue)	Biosima Model	Volume(µL)	Sterile	Low-Retention	Filtering	Packing
 Tip Marks: 200µL, 500µL	T-1000B-B	100-1000				1000 pcs/pack, 5 packs/carton
	T-1000B-BL	100-1000		✓		1000 pcs/pack, 5 packs/carton
	T-1000B-R	100-1000				96 pcs/rack,10 racks /box, 5 boxes/carton
	T-1000B-RL	100-1000		✓		96 pcs/rack,10 racks /box, 5 boxes/carton
	T-1000B-RS	100-1000	✓			96 pcs/rack,10 racks /box, 5 boxes/carton
	T-1000B-RSL	100-1000	✓	✓		96 pcs/rack,10 racks /box, 5 boxes/carton
	T-1000B-E	100-1000				96 pcs/rack, 5 racks/box, 10 boxes/carton
	T-1000B-EL	100-1000		✓		96 pcs/rack, 5 racks/box, 10 boxes/carton
	T-1000B-F	100-1000			✓	96 pcs/rack, 5 racks/box, 10 boxes/carton

1500µL Tips	Biosima Model	Volume(µL)	Sterile	Low-Retention	Filtering	Packing
 Tip Marks: 500µL, 1000µL, 1500µL	TF-1500-B	100-1500			✓	1000 pcs/pack, 5 packs/carton
	TF-1500-BL	100-1500		✓	✓	1000 pcs/pack, 5 packs/carton
	TF1500-R	100-1500			✓	96 pcs/rack,10 racks /box, 5 boxes/carton
	TF-1500-RL	100-1500		✓	✓	96 pcs/rack,10 racks /box, 5 boxes/carton
	TF-1500-RS	100-1500	✓		✓	96 pcs/rack,10 racks /box, 5 boxes/carton
	TF-1500-RSL	100-1500	✓	✓	✓	96 pcs/rack,10 racks /box, 5 boxes/carton

Compatible with a variety of pipettors, suitable for pipetting micro reagent from deep vessels.

Integra Pipette Tips

Integra Tips	Biosima Model	Volume(µL)	Sterile	Low-Retention	Filtering	Packing
 Tip Marks: 12.5µL	IG-384-12.5-RS	12.5	✓			384 pcs/rack, 5 racks/box, 4 boxes/carton
	IG-384-12.5-HS	12.5	✓			384 pcs/rack, 5 racks/box, 4 boxes/carton
 Tip Marks: 12.5µL	IG-384-12.5L-RS	12.5	✓			384 pcs/rack, 5 racks/box, 4 boxes/carton
	IG-384-12.5L-HS	12.5	✓			384 pcs/rack, 5 racks/box, 4 boxes/carton
 Tip Marks: 125µL	IG-384-125-RS	125	✓			384 pcs/rack, 5 racks/box, 4 boxes/carton
 Tip Marks: 300µL	IG-96-300-RS	300	✓			96 pcs/rack, 10 racks/box, 5 boxes/carton
	IG-96-300-RSL	300	✓	✓		96 pcs/rack, 10 racks/box, 5 boxes/carton
	IGF-96-300-RS	300	✓		✓	96 pcs/rack, 10 racks/box, 5 boxes/carton
	IGF-96-300-RSL	300	✓	✓	✓	96 pcs/rack, 10 racks/box, 5 boxes/carton
 Tip Marks: 1250µL	IG-96-1250-RS	1250	✓			96 pcs/rack, 10 racks/box, 5 boxes/carton
	IG-96-1250-RSL	1250	✓	✓		96 pcs/rack, 10 racks/box, 5 boxes/carton
	IGF-96-1250-RS	1250	✓		✓	96 pcs/rack, 10 racks/box, 5 boxes/carton
	IGF-96-1250-RSL	1250	✓	✓	✓	96 pcs/rack, 10 racks/box, 5 boxes/carton

Compatible with Integra pipettors, suitable for pipette micro liquid from deep vessels.

For Rainin Pipettor

1



20µL Tips

2



200µL Tips

3





300µL Tips

4





1000µL Tips



20µL Tips	Biosima Model	Volume(µL)	Sterile	Low-Retention	Filtrating	Packing
 φ6.2mm φ0.7mm 49.4mm Tip Marks: 2µL, 10µL, 20µL	LTS-20-B	0.5-20				1000 pcs/pack, 10 packs/carton
	LTS-20-BL	0.5-20		✓		1000 pcs/pack, 10 packs/carton
	LTS-20-R	0.5-20				96 pcs/rack, 10 racks /box, 5 boxes/carton
	LTS-20-RL	0.5-20		✓		96 pcs/rack, 10 racks /box, 5 boxes/carton
	LTS-20-RS	0.5-20	✓			96 pcs/rack, 10 racks /box, 5 boxes/carton
	LTS-20-RSL	0.5-20	✓	✓		96 pcs/rack, 10 racks /box, 5 boxes/carton
	LTS-20-E	0.5-20				96 pcs/rack, 10 racks /box, 10 boxes/carton
	LTS-20-EL	0.5-20		✓		96 pcs/rack, 10 racks /box, 10 boxes/carton
	LTSF-20-B	0.5-20			✓	1000 pcs/pack, 10 packs/carton
	LTSF-20-BL	0.5-20		✓	✓	1000 pcs/pack, 10 packs/carton
	LTSF-20-RS	0.5-20	✓		✓	96 pcs/rack, 10 racks /box, 5 boxes/carton
	LTSF-20-RSL	0.5-20	✓	✓	✓	96 pcs/rack, 10 racks /box, 5 boxes/carton
200µL Tips	Biosima Model	Volume(µL)	Sterile	Low-Retention	Filtrating	Packing
 φ7.3mm φ0.8mm 49.9mm Tip Marks: 10µL, 50µL, 100µL, 200µL	LTS-200-B	10-200				1000 pcs/pack, 10 packs/carton
	LTS-200-BL	10-200		✓		1000 pcs/pack, 10 packs/carton
	LTS-200-R	10-200				96 pcs/rack, 10 racks /box, 5 boxes/carton
	LTS-200-RL	10-200		✓		96 pcs/rack, 10 racks /box, 5 boxes/carton
	LTS-200-RS	10-200	✓			96 pcs/rack, 10 racks /box, 5 boxes/carton
	LTS-200-RSL	10-200	✓	✓		96 pcs/rack, 10 racks /box, 5 boxes/carton
	LTS-200-E	10-200				96 pcs /layer, 10 layers /box, 10 boxes/carton
	LTS-200-EL	10-200		✓		96 pcs /layer, 10 layers /box, 10 boxes/carton
	LTSF-200-B	10-200			✓	1000 pcs/pack, 10 packs/carton
	LTSF-200-BL	10-200		✓	✓	1000 pcs/pack, 10 packs/carton
	LTSF-200-RS	10-200	✓		✓	96 pcs/rack, 10 racks /box, 5 boxes/carton
	LTSF-200-RSL	10-200	✓	✓	✓	96 pcs/rack, 10 racks /box, 5 boxes/carton

Compatible with Rainin, For LiteTouch System, XLS pipette.

30µL Tips	Biosima Model	Volume(µL)	Sterile	Low-Retention	Filtering	Packing
	LTS-300-B	10-300				1000 pcs /pack, 10 packs /carton
	LTS-300-BL	10-300		✓		1000 pcs /pack, 10 packs /carton
	LTS-300-R	10-300				96 pcs /rack,10 racks /box, 5 boxes /carton
	LTS-300-RL	10-300		✓		96 pcs /rack,10 racks /box, 5 boxes /carton
	LTS-300-RS	10-300	✓			96 pcs /rack,10 racks /box, 5 boxes /carton
	LTS-300-RSL	10-300	✓	✓		96 pcs /rack,10 racks /box, 5 boxes /carton
	LTS-300-E	10-300				96 pcs /layer, 8 layers /box, 10 boxes /carton
	LTS-300-EL	10-300		✓		96 pcs /layer, 8 layers /box, 10 boxes /carton
	LTSF-300-B	10-300			✓	1000 pcs /pack, 10 packs /carton
	LTSF-300-BL	10-300		✓	✓	1000 pcs /pack, 10 packs /carton
	LTSF-300-RS	10-300	✓		✓	96 pcs /rack,10 racks /box, 5 boxes /carton
	LTSF-300-RSL	10-300	✓	✓	✓	96 pcs /rack,10 racks /box, 5 boxes /carton
Tip Marks: 10µL, 50µL, 100µL						
1000µL Tips	Biosima Model	Volume(µL)	Sterile	Low-Retention	Filtering	Packing
	LTS-1000-B	100-1000				1000 pcs /pack, 10 packs /carton
	LTS-1000-BL	100-1000		✓		1000 pcs /pack, 10 packs /carton
	LTS-1000-R	100-1000				96 pcs /rack,10 racks /box, 5 boxes /carton
	LTS-1000-RL	100-1000		✓		96 pcs /rack,10 racks /box, 5 boxes /carton
	LTS-1000-RS	100-1000	✓			96 pcs /rack,10 racks /box, 5 boxes /carton
	LTS-1000-RSL	100-1000	✓	✓		96 pcs /rack,10 racks /box, 5 boxes /carton
	LTS-1000-E	100-1000				96 pcs /layer, 8 layers /box, 10 boxes /carton
	LTS-1000-EL	100-1000		✓		96 pcs /layer, 8 layers /box, 10 boxes /carton
	LTSF-1000-B	100-1000			✓	1000 pcs /pack, 5 packs /carton
	LTSF-1000-BL	100-1000		✓	✓	1000 pcs /pack, 5 packs /carton
	LTSF-1000-RS	100-1000	✓		✓	96 pcs /rack,10 racks /box, 5 boxes /carton
	LTSF-1000-RSL	100-1000	✓	✓	✓	96 pcs /rack,10 racks /box, 5 boxes /carton
Tip Marks: 100µL						

Compatible with Rainin, For LiteTouch System, XLS pipette.

Pipette Tip Package



Rack Package

- ✓ **One-Handed Performance**
To improve the working efficiency.
- ✓ **Translucent Lid**
Different color for different volume, easy to recognize.
- ✓ **Individual Rack Package**
individual package to prevent pollution.

Reloading Package

- ✓ **To Save Space**
There are 8 or 10 layers tips in every package , and it can be reloaded.
- ✓ **Easy To Operate**
It can reload the individual rack in some seconds by the reloading package.
- ✓ **Environmental Friendly**
The package way can save 60% plastic usage.



Individually Rack
Easy to use
96 pcs /rack



Stacked Box Packaging
Save space and reduce plastic waste by 60%
96 pcs /layer



Bagged Packing
Save space and reduce plastic waste by 85%
1000 pcs /bag

Variety of sizes, types and packaging are available to meet needed.

■ 20 µL ■ 200 µL ■ 300 µL ■ 1000 µL

A vertical pipette is shown dispensing a blue liquid into a multi-well plate. The pipette has a grey body and a blue tip. The liquid is being dispensed into one of the wells of the plate. The text "Tips Adaptation Form" is overlaid at the top of the image.


DLAB

$\sqrt{r}$ r =recommend

✓r r=recommend

CV < 3.5%

Biosima tips ensures consistency and repeatability of results between operators.



Deep Well Plate



SBS International Standard

All products conform to SBS international standard, are compatible with most of pipette workstations.



Imported Virgin Material

Imported medical grade polypropylene (PP).



High Temperature Sterilization

Autoclavable at 121°C/15psi.



Planeness

Smooth surface and uniform spacing between wells.



Non - Pyrogenic

Non - pyrogen, DNase/RNase Free.



Chemical Resistance

Resistant to chemical reagent during use and storage.



Traceable

All product package is coded for traceability.

The Applications Of Deep Well Plate

> Storage Samples


It is application to store samples, perfectly replace the 1.5ml centrifuge tube, with large storage capacity, saving more space, working well in -80 °C.

> Sample treatment

It achieves high-throughput operation of biological samples with pipette and high-throughput automatic liquid operation instruments and software.

> Operation

It is applicable to all kinds of automatic equipment, can be directly put in the sample bin and make operation more efficiency.



DEEP WELL PLATE

Bottom Shape Of Deep Well Plate



U Bottom

Better mixing of samples.



Tapered Bottom / V Bottom

Better the sample recovery and less liquid residue. Increased heating wrap area to make heating more efficient.




V Bottom


Almost no Invalid volume in the tube during liquid treatment.



U Bottom



Tapered Bottom



V Bottom




U Bottom



V Bottom







Flat Bottom




V-Tapered Bottom

PCR Plate

			
PCR-0208	PCR-96-HN-C	PCR-96-HS-C	PCR-96R-HS-C
PCR Tube Strips	96 PCR Plate (No Profile)	96 PCR Plate (Half Profile)	96 PCR Plate (Half Profile, Removable)
8 Holes (0.2 mL)	96 Holes (0.2 mL)	96 Holes (0.2 mL)	96 Holes (0.2 mL)
			
71X8X20.6 mm	119.8X79.7X21 mm	126.3X86X21.1 mm	127.76X85.48X20.7 mm
125 pcs/bag, 10 bags/carton	10 pcs/bag, 5 bags/carton	10 pcs/bag, 5 bags/carton	10 pcs/bag, 5 bags/carton


A full range of PCR plates, printing services can be provided.






PCR Cover (Sealing film/Silicone Mat) :









	
F-101	M-PCR-RD-96
Adaptation, full series of PCR plates	Adaptation, full series of PCR plates
Features: PCR fluorescence quantification	Features: Round hole, pierceable
130X80 mm	113.5X77.5 mm
125 pcs/bag, 10 bags/carton	10 pcs/bag, 5 bags/carton

DEEP WELL PLATE

Round Well Plate











			
P-0.36-96V	P-0.4-96F	P-0.4-96U	P-0.45-RDV-96
V Bottom	Flat Bottom	U Bottom	V Bottom
96 Holes	96 Holes	96 Holes	96 Holes
0.36 mL	0.4 mL	0.4 mL	0.45 mL
			
127.8X85.5X14.8 mm	127.8X85.5X14.8 mm	127.8X85.5X14.8 mm	127.15X85.2X14.5 mm
10 pcs/bag, 100 pcs/carton	10 pcs/bag, 100 pcs/carton	10 pcs/bag, 100 pcs/carton	10 pcs/bag, 100 pcs/carton

			
P-0.5-RD-96	P-1.0-RD-96	P-1.0-RDV-96	P-1.0-RDB-96
U Bottom	U Bottom, Independent Well	V Bottom	U Bottom
96 Holes	96 Holes	96 Holes	96 Holes
0.5 mL	1.0 mL	1.0 mL	1.0 mL
			
127.15X85.2X14.5 mm	127.8X85.5X33 mm	127.8X85.5X33 mm	127.76X85.48X41.7 mm
10 pcs/bag, 100 pcs/carton	5 pcs/bag, 50 pcs/carton	5 pcs/bag, 50 pcs/carton	5 pcs/bag, 50 pcs/carton







DEEP WELL PLATE








Round Well Plate

			
P-1.1-RD-96	P-1.1-RDV-96	P-2.0-RDA-96	P-2.0-RDV-96
U Bottom, Shared Well	V Bottom, Shared Well	U Bottom	V Bottom
96 Holes	96 Holes	96 Holes	96 Holes
1.1 mL	1.1 mL	2.0 mL	2.0 mL
			
127.5X85.35X27 mm	127.5X85.35X27 mm	127.5X85.35X44.3 mm	127.5X85.35X44.3 mm
5 pcs/bag, 50 pcs/carton	5 pcs/bag, 50 pcs/carton	5 pcs/bag, 50 pcs/carton	5 pcs/bag, 50 pcs/carton

		
P-3.5-RD-48	P-3.5-RDA-48	P-C
U Bottom, Irregular Arrangement	U Bottom, Regular Arrangement	Transparent Plate Cover
48 Holes	48 Holes	/
3.5 mL	3.5 mL	Adaptation, full series of PCR plates
		/
127.8X85.5X43.6 mm	127.5X85.2X51.1 mm	129.2X86X9.3 mm
5 pcs/bag, 50 pcs/carton	5 pcs/bag, 50 pcs/carton	20 pcs/bag, 200 pcs/carton






Square Well Plate






			
P-0.5-SQV-96	P-1.0-SQ-96	P-1.0-SQV-96	P-1.6-SQ-96
Tapered Bottom	U Bottom	Tapered Bottom	U Bottom
96 Holes	96 Holes	96 Holes	96 Holes
0.5 mL	1.0 mL	1.0 mL	1.6 mL
			
127.5X85.4X15.1 mm	127.8X85.5X24.4 mm	127.4X85.5X24.4 mm	127.5X85.4X31.6 mm
10 pcs/bag, 100 pcs/carton	10 pcs/bag, 100 pcs/carton	10 pcs/bag, 100 pcs/carton	5 pcs/bag, 50 pcs/carton

			
P-1.6-SQV-96	P-2.2-SQ-96	P-2.2-SQG-96	P-2.2-SQV-96
V Bottom	U Bottom	U Bottom, Profile Concave	Tapered Bottom, Frosted Shell
96 Holes	96 Holes	96 Holes	96 Holes
1.6 mL	2.2 mL	2.2 mL	2.2 mL
			
127.5X85.4X31.6 mm	127.5X85.4X41.6 mm	127X85X42.5 mm	127.4X85.5X44.2 mm
5 pcs/bag, 50 pcs/carton	5 pcs/bag, 50 pcs/carton	5 pcs/bag, 50 pcs/carton	5 pcs/bag, 50 pcs/carton


DEEP WELL PLATE

Square Well Plate





			
P-2.2-SQVT-96	P-2.3-SQV-96	P-2.3-SQC-96	P-10-SQV-24
Tapered Bottom, Transparent Shell	V Bottom	Diamond Cone Bottom	Tapered Bottom
96 Holes	96 Holes	96 Holes	24 Holes
2.2 mL	2.3 mL	2.3 mL	10 mL
			
127.4X85.5X44.2 mm	127.5X85.4X41.7 mm	127.3X85.2X44.3 mm	127.5X85.5X44 mm
5 pcs/bag, 50 pcs/carton	5 pcs/bag, 50 pcs/carton	5 pcs/bag, 50 pcs/carton	5 pcs/bag, 50 pcs/carton

			
P-15-SQV-24	P-4.6-SQ-48	P-120-SQ-384	P-240-SQ-384
Tapered Bottom	U Bottom	V Bottom	V Bottom
24 Holes	48 Holes	384 Holes	384 Holes
15 mL	4.6 mL	120 µL	240 µL
			
127.2X85.3X63.2 mm	127.8X85.5X41.6 mm	127.5X85.35X14.6 mm	127.5X83.35X24.9 mm
5 pcs/bag, 50 pcs/carton	5 pcs/bag, 50 pcs/carton	10 pcs/white box, 100 pcs/carton	10 pcs/white box, 50 pcs/carton





Silicone Mat For Well Plate





DEEP WELL PLATE


			
M-RD-96-U	M-2mL-RD	M-2mL-RD-C	M-RD-3.5mL
Adaptation 1.0 mL round well plate	Adaptation 2.0 mL round well plate	Adaptation 2.0 mL round well plate	Adaptation 3.5 mL round well plate
Features Round hole, can be pierced	Features Round hole, can be pierced	Features Cross cut, pierceable	Features Round hole, pierceable
10 pcs/bag, 5 bags/carton	10 pcs/bag, 5 bags/carton	10 pcs/bag, 5 bags/carton	10 pcs/bag, 5 bags/carton

Silicone Mat

			
M-SQ-96	M-SQ-96-P	M-SQ-96-C	M-SQ-24
Adaptation 96 square well plate	Adaptation 96 square well plate	Adaptation 96 square well plate	Adaptation 24 square well plate
Features Square hole	Features Square hole, pierceable	Features Square hole, cross cut	Features Square hole
10 pcs/bag, 5 bags/carton	10 pcs/bag, 5 bags/carton	10 pcs/bag, 5 bags/carton	10 pcs/bag, 5 bags/carton

	
M-SQ-48	M-SQ-384
Adaptation 48 square well plate	Adaptation 384 square well plate
Features Square hole	Features Square hole
10 pcs/bag, 50 pcs/carton	10 pcs/bag, 50 pcs/carton





Sealing Film








DEEP WELL PLATE




Biosima Model	Dimensions(mm)	Description	Packing
F-001	144x87	Adhesive sealing film, Suit for Deep well plates and PCR	500 pcs/bag, 20 bags /case
F-002	125x82	Embossed, Suit for deep Well plate, PCR	600 pcs/bag, 12 bags /case
F-004	127x78	Imported, Pierceable, Glossy, For deep Well plate, PCR	100 pcs/bag, 5 bags /case
F-005	125x83	Glossy, Thick, Golden, For deep Well plate, PCR	600 pcs/bag, 16 bags /case
F-007	124x78	Glossy, Thick, Golden, For deep Well plate, PCR	100 pcs/bag, 5 bags /case
F-009	124x78	Imported, Glossy, Thick, For deep Well plate, PCR	100 pcs/bag, 5 bags /case
F-012	125x83	Glossy, Thick, For deep Well plate, PCR	100 pcs/bag, 5 bags /case
F-101	141x77	Fluorescent quantitation sealing film	100 pcs/bag, 5 bags /case
F-103	140x80	Piercable, PCR-Sensitive Membrane	100 pcs/bag, 5 bags /case
F-104	100x55 (48) 100x25 (24)	PCR-Sensitive Membrane For 24 Well/48 Well plate	100 pcs/bag, 5 bags /case
F-105	140x83 (Back Paper) 137x80 (Top Paper)	High transmittance adhesive sealing film	100 pcs/bag, 5 bags /case
F-106	140x83 (Back Paper) 133x80 (Top Paper)	High transmittance adhesive sealing film	100 pcs/bag, 5 bags /case
F-107	140x83 (Back Paper) 137x80 (Top Paper)	Resistant to 150°C high transmittance sealing film	100 pcs/bag, 5 bags /case
F-108	128x90 (Full Plate) 90x9.2 (Single Strip)	detachable sealing film for ELISA Plate	100 pcs/bag, 5 bags /case
F-109	140x83 (Back Paper) 137x80 (Top Paper)	High transmittance sealing film for culture plate	100 pcs/bag, 5 bags /case
F-110	140x83 (Back Paper) 137x80 (Top Paper)	High transmittance sealing film for culture plate	100 pcs/bag, 5 bags /case
F-111	140x83 (Back Paper) 137x80 (Top Paper)	High transmittance sealing film for culture plate	100 pcs/bag, 5 bags /case

Tip Comb And Pc-Cover

			
P-MT-8U	J-MT-8V	D-MT-8U	SS-MT-8U
U Bottom	V Bottom, Snap	U Bottom, Snap	U Bottom, Snap
8 Holes	8 Holes	8 Holes	8 Holes
Adaptation 2.2 mL square well plate	Adaptation 2.2 mL square well plate	Adaptation 2.2 mL square well plate	Adaptation 2.2 mL square well plate
			
11.45X44.7 mm	12X45.5 mm	12X46.1 mm	11.8X45.8 mm
200 sets/carton	200 sets/carton	200 sets/carton	200 sets/carton

			
P-MT-24	P-MT-24L	P-MT-96V	P-C
V Bottom	V Bottom	V Bottom	Transparent Plate Cover
24 Holes	24 Holes	96 Holes	/
Adaptation 10 mL square well plate	Adaptation 15 mL square well plate	Adaptation 2.2 mL square well plate	Adaptation full range of well plates
			/
126.5X84.9X44.2 mm	127.4X85.3X62.2 mm	127.3X85.3X42.3 mm	129.2X86X9.3 mm
50 sets/carton	50 sets/carton	50 sets/carton	20 pcs/bag, 200 pieces/carton

Compatible With TANBEAD

		
	P-SQ-84R	P-MT-12U
	Square Well Pate, Detachable	Snap, PP
	84 Holes	12 Holes
	Suit for TANBEAD pipetting workstation	Suit for TANBEAD square well pate
		
	152.3X81.8X52.7 mm	12X56.2 mm
	16 pcs/bag, 30 bags/carton	15 pcs/bag, 4 bags/carton



Suit for TANBEAD SLA-DI4800
pipetting workstation.

DEEP WELL PLATE

Reservoir

The Application Of Reservoir



The reservoir has standard SBS size, which is suitable for Plate reactor and processing system, and can process reagents in a simple and economical way on the automatic system. The reservoir is disposable, the built-in groove at the bottom of the reservoir can reduce the waste of reagent and ensure the consistency of washing liquid.

RESERVIOR



SBS International Standard

All products conform to SBS international standard, are compatible with most of pipette workstations.



Imported Virgin Material

Imported medical grade polypropylene (PP).



High Temperature Sterilization

Autoclavable at 121°C/15psi.



Non - Pyrogenic

Non - pyrogen, DNase/RNase Free.



Chemical Resistance

Resistant to chemical reagent during use and storage.



Low Heavy Metal Content

High quality and low heavy metal content.



Traceable

All product package is coded for traceability.


Smooth Pipetting

In the process of automatic liquid transfer, the liquid cannot be taken evenly due to the surface tension in the flat bottom container. But reservoir with groove can overcome the tension and reduce the liquid residue.



Flat Bottom Container









It's hard to transfer liquid evenly due to the liquid tension.










Reservoir





The liquid is evenly distributed in the reservoir.






Reservoir

			
RES-12-8MW	RES-32-8MW	RES-180-8MW	RES-8-12MW
Multi Well, Low Profile	Multi Well, High Profile	Multi Well, Middle Profile	Multi Well, Low Profile
8 Channel	8 Channel	8 Channel	12 Channel
12 mL	32 mL	22 mL	8 mL
			
127.76X85.48X19.2 mm	127.76X85.48X44.5 mm	127.76X85.48X31.5 mm	127.76X85.48X19.2 mm
10 pcs/bag, 5 bags/carton	5 pcs/bag, 50 bags/carton	5 pcs/bag, 50 bags/carton	10 pcs/bag, 5 bags/carton

			
RES-22-12MW	RES-180-12MW	RES-100-F	RES-190-F
Multi Well, High Profile	Multi Well, Middle Profile	Single Well, Low Profile	Single Well, Middle Profile
12 Channel	12 Channel	Flat Channel	Flat Channel
22 mL	15 mL	100 mL	190 mL
			
127.76X85.48X44.5 mm	127.76X85.48X31.5 mm	127.76X85.48X19.2 mm	127.76X85.48X31.5 mm
5 pcs/bag, 50 bags/carton	5 pcs/bag, 50 bags/carton	10 pcs/bag, 5 bags/carton	5 pcs/bag, 50 bags/carton

Reservoir

			
RES-100-96SW	RES-195-96SW	RES-290-8SW	RES-290-12SW
Single Well, Low Profile	Single Well, Middle Profile	Single Well, High Profile	Single Well, High Profile
96 Channel	96 Channel	8 Channel	12 Channel
100 mL	195 mL	290 mL	290 mL
			
127.76X85.48X19.2 mm	127.76X85.48X31.5 mm	127.76X85.48X44.5 mm	127.76X85.48X44.5 mm
10 pcs/bag, 5 bags/carton	5 pcs/bag, 50 bags/carton	5 pcs/bag, 50 bags/carton	5 pcs/bag, 50 bags/carton

			
RES-290-96SW	RES-100-384SW	RES-185-384SW	RES-290-384SW
Single Well, High Profile	Single Well, Low Profile	Single Well, Middle Profile	Single Well, High Profile
96 Channel	384 Channel	384 Channel	384 Channel
290 mL	100 mL	185 mL	290 mL
			
127.76X85.48X44.5 mm	127.76X85.48X19.2 mm	127.76X85.48X31.5 mm	127.76X85.48X44.5 mm
5 pcs/bag, 50 bags/carton	10 pcs/bag, 5 bags/carton	5 pcs/bag, 50 bags/carton	5 pcs/bag, 50 bags/carton

Reagent Tank



10 mL Reagent Slot	Biosima Model	Volume(mL)	Type	Packing
	RES-PS10	10	Bulk	300 pc/bag
	RES-PS10-1S	10	Individual Pack, Sterilization	1 pc/bag, 100 bags/case
	RES-PS10-5S	10	5 pc/bag, Sterilization	5 pc/bag, 40 bags/case


25 mL Reagent Slot	Biosima Model	Volume(mL)	Type	Packing
	RES-PS25	25	Bulk	100 pc/bag
	RES-PS25-1S	25	Individual Pack, Sterilization	1 pc/bag, 100 bags/case
	RES-PS25-5S	25	5 pc/bag, Sterilization	5 pc/bag, 40 bags/case
50 mL Reagent Slot	Biosima Model	Volume(mL)	Type	Packing
	RES-PS50	50	Bulk	100 pc/bag
	RES-PS50-1S	50	Individual Pack, Sterilization	1 pc/bag, 100 bags/case
	RES-PS50-5S	50	5 pc/bag, Sterilization	5 pc/bag, 40 bags/case
65 mL Reagent Slot	Biosima Model	Volume(mL)	Type	Packing
	RES-PS65	65	Bulk	100 pc/bag
	RES-PS65-1S	65	Individual Pack, Sterilization	1 pc/bag, 100 bags/case
	RES-PS65-5S	65	5 pc/bag, Sterilization	5 pc/bag, 40 bags/case
100 mL Reagent Slot	Biosima Model	Volume(mL)	Type	Packing
	RES-PS100	100	Bulk	100 pc/bag
	RES-PS100-1S	100	Individual Pack, Sterilization	1 pc/bag, 100 bags/case
	RES-PS100-5S	100	5 pc/bag, Sterilization	5 pc/bag, 40 bags/case

RESERVIROR

Tube & Column

The Application Of Centrifuge Tube

Column nucleic acid extraction is composed of adsorption column and collection tube. This method is suitable for extracting plasmids, viruses, DNA, RNA and so on.



Imported Virgin Material

Imported medical grade polypropylene (PP).



High Temperature Sterilization

Autoclavable at 121°C/15psi.



Refrigerated

The minimum enduring temperature can reach - 80 °C.



Centrifugal Force

Can withstand a maximum centrifugal force of 22000xg without damage or deformation.



Low Residue

Less residual liquid, less wall hanging phenomenon.




Chemical Resistance

Resistant to chemical reagent during use and storage.



Traceable

All product package is coded for traceability.



1 Black Scale
Easy to read data.

2 Convenience Mark
Tube body with large marking area.

3 Precise Size
Compatible with all brands of centrifuge rotors and mixers.












4 High Quality Materials
Corrosion resistance, high centrifugal stability.

5 Sterile And Pollution-free
DNase-free, RNase-free, pyrogen-free.



Tube/Column




			
T-001	T-002	T-003	T-004
Centrifuge Tube with long cap	Centrifuge Tube with long cap	Spin Preparation Column with ring	Spin Gel Column with ring
1.5 mL	2.0 mL	0.7 mL	0.7 mL
Centrifugal Force 22000 xg	Centrifugal Force 22000 xg	Centrifugal Force 22000 xg	Centrifugal Force 22000 xg
50 pcs/bag, 100 bags/carton	50 pcs/bag, 100 bags/carton	1000 pcs/bag, 10 bags/carton	1000 pcs/bag, 10 bags/carton
			
T-005	T-007	T-008	T-009
micro centrifuge tube with cap	Purification Spin Column Flat bottom/with cover and ring	micro centrifuge tube no cap, non-graduated	Purification Spin Column Including pressure ring
1.5 mL	0.7 mL	2.0 mL	0.7 mL
Centrifugal Force 22000 xg	Centrifugal Force 22000 xg	Centrifugal Force 22000 xg	Centrifugal Force 22000 xg
500 pcs/bag, 10 bags/carton	1000 pcs/bag, 10 bags/carton	1000 pcs/bag, 15 bags/carton	1000 pcs/bag, 10 bags/carton
			
CT-1.5	CT-1.5-L	CT-2.0	CT-2.0-L
Centrifuge Tube Double explosion proof buckle	Double explosion proof buckle Low retention	Centrifuge Tube Double explosion proof buckle	Double explosion proof buckle Low retention
1.5 mL	1.5 mL	2.0 mL	2.0 mL
Centrifugal Force 22000 xg	Centrifugal Force 22000 xg	Centrifugal Force 22000 xg	Centrifugal Force 22000 xg
500 pcs/bag, 10 bags/carton	500 pcs/bag, 10 bags/carton	500 pcs/bag, 10 bags/carton	500 pcs/bag, 10 bags/carton





Tube/Column

			
CT-15-S	CT-15-RS	CT-15SS-S	CT-15SS-RS
Centrifuge Tube sterilized	Centrifuge Tubes Rack mounted, sterilized	Centrifuge Tubes Can stand, sterilized	Stand mounted Can stand, sterilized
15 mL	15 mL	15 mL	15 mL
Centrifugal Force 12000 xg	Centrifugal Force 12000 xg	Centrifugal Force 12000 xg	Centrifugal Force 12000 xg
25 pcs/bag, 20 bags/carton	25 pcs/bag, 20 bags/carton	25 pcs/bag, 20 bags/carton	25 pcs/bag, 20 bags/carton

			
CT-25-S	CT-25SS-S	CT-50-S	CT-50-RS
Centrifuge Tube sterilized	Centrifuge Tubes Can stand, sterilized	Centrifuge Tube sterilized	Centrifuge Tubes Rack mounted, sterilized
25 mL	25 mL	50 mL	50 mL
Centrifugal Force 12000 xg	Centrifugal Force 12000 xg	Centrifugal Force 12000 xg	Centrifugal Force 12000 xg
25 pcs/bag, 8 bags/carton	25 pcs/bag, 8 bags/carton	25 pcs/bag, 20 bags/carton	25/rack, 12/box

			
CT-50SS-S	CT-50SS-RS	CT-50SSH-S	CT-50SSH-RS
Centrifuge Tubes Can stand, sterilized	Stand mounted Can stand, sterilized	Can stand, hexagonal bottom sterilization	Bracket mounted, can stand Hexagonal bottom, sterilized
50 mL	50 mL	50 mL	50 mL
Centrifugal Force 12000 xg	Centrifugal Force 12000 xg	Centrifugal Force 12000 xg	Centrifugal Force 12000 xg
25 pcs/bag, 20 bags/carton	25/rack, 12/box	25/rack, 12/box	25/rack, 12/box

Tube/Column Box


			
CTR-15	CTR-50	DWP-24PPK-C	DWP-24PPK-C-S
Centrifuge Tube Holder	Centrifuge Tube Holder	Centrifuge Tube Box	Centrifuge Tube Box
Adaptation 15 mL centrifuge tube /	Adaptation 50 mL centrifuge tube /	Adaptation T-005,CT-1.5,CT-2.0 /	Adaptation T-005,CT-1.5,CT-2.0 /
/	/	5 pcs/bag, 50 pcs/carton	5 pcs/bag, 50 pcs/carton

Tube/Column Package

		
Individually Rack	Bracket Packaging	Bagged Packaging
Easy to use	Save space and reduce plastic waste by 60%	
24/box	25/rack	200 pcs/bag, 500 pcs/bag

Variety of sizes, types and packaging are available to meet needed.

Reagent Bottle



Type	Diameter (mm)				Height (mm)				Bottle Cap (mm)			
	A	Deviate	B	Deviate	C	Deviate	D	Deviate	E	Deviate	F	Deviate
B-8mL	ø24.5	±0.5	ø17.5	±0.5	42.6	±1.0	44	±1.0	ø24.5	±0.5	15	±0.5
B-15mL	ø24.5	±0.5	ø17.5	±0.5	55.2	±1.0	57.5	±1.0	ø24.5	±0.5	15	±0.5
B-30mL	ø33.5	±0.5	ø25	±0.5	58.7	±1.0	61.5	±1.5	ø32.0	±0.5	15.5	±0.5
B-60mL	ø39	±1.0	ø25	±0.5	80	±1.0	83	±1.5	ø32.0	±0.5	15.5	±0.5
B-125mL	ø50	±1.0	ø32	±0.5	94	±1.0	99	±1.5	ø42	±0.5	17	±0.8
B-250mL	ø60.5	±1.0	ø36.5	±0.5	126	±1.0	130	±1.5	ø47	±0.5	16	±0.8
B-500mL	ø71.5	±1.0	ø47	±0.5	161	±2.0	166	±2.0	ø57	±1.0	17	±1.0
B-1000mL	ø91	±2.0	ø56.5	±0.5	192	±2.0	199	±2.0	ø68	±1.5	18	±1.5
B-1000mL-N	ø92	±2.0	ø34	±0.5	211	±2.0	215	±2.0	ø42	±1.5	30	±1.5

Reagent Bottle



B-8/B-8-B

Clear/Brown

8 mL

1000 sets/box



B-15/B-15-B

Clear/Brown

15 mL

1000 sets/box



B-30/B-30-B

Clear/Brown

30 mL

400 sets/box



B-60/B-60-B

Clear/Brown

60 mL

650 sets/box



B-125/B-125-B

Clear/Brown

125 mL

350 sets/box



B-250/B-250-B

Clear/Brown

250 mL

200 sets/box



B-500/B-500-B

Clear/Brown

500 mL

100 sets/box



B-1000/B-1000-B

Clear/Brown

1000 mL

50 sets/box



B-1000-N/B-1000-B-N

Clear/Brown

1000 mL

50 sets/box

Serological Pipettes



USP Raw Material
USP grade VI polystyrene (PS) raw material.



High Temperature Sterilization
Autoclavable at 121°C/15psi.



Color Ring For Easier Identification
Different colors correspond to different specification, easy to distinguish.



Bidirectional Scale
There are ascending and descending scales to meet the requirements of simultaneous sampling and sampling reduction.



Non - Pyrogenic
Non – pyrogen, DNase/RNase Free.



Chemical Resistance
Resistant to chemical reagent during use and storage.



Low Heavy Metal Content
High quality and low heavy metal content.



Traceable
All product package is coded for traceability.

Serological Pipettes



The product is made of polystyrene (PS) conforming to USP grade VI. It can fit most of the pipettors on the market. This series of products are mainly used to accurately measure or transfer a certain volume of liquid.

Sterilization Packaging



Paper-Plastic Bags
Individually packaged in paper-plastic bags.




Plastic Bags
Individually packaged in plastic bags.

1mL Serological Pipettes (Sterilization)




Biosima Model	Volume	Scale	Color Ring	Packing
SP01001	1 mL	1/100 mL		Paper-Plastic Bags, 50 pcs/bag, 20 bags/case
SP02001	1 mL	1/100 mL		Plastic Bags, 50 pcs/bag, 20 bags/case

2mL Serological Pipettes (Sterilization)



Biosima Model	Volume	Scale	Color Ring	Packing
SP01002	2 mL	1/50 mL		Paper-Plastic Bags, 50 pcs/bag, 20 bags/case
SP02002	2 mL	1/50 mL		Plastic Bags, 50 pcs/bag, 20 bags/case

5mL Serological Pipettes (Sterilization)



Biosima Model	Volume	Scale	Color Ring	Packing
SP01005	5 mL	1/10 mL		Paper-Plastic Bags, 50 pcs/bag, 10 bags/case
SP02005	5 mL	1/10 mL		Plastic Bags, 50 pcs/bag, 10 bags/case


10mL Serological Pipettes (Sterilization)



Biosima Model	Volume	Scale	Color Ring	Packing
SP01010	10 mL	1/10 mL		Paper-Plastic Bags, 50 pcs/bag, 10 bags/case
SP02010	10 mL	1/10 mL		Plastic Bags, 50 pcs/bag, 10 bags/case


Compatible with all kinds of pipettes, suitable for measuring a certain volume of liquid and vacuum suction.

25mL Serological Pipettes (Sterilization)




Biosima Model	Volume	Scale	Color Ring	Packing
SP01025	25 mL	2/10 mL		Paper-Plastic Bags, 25 pcs/bag, 10 bags/case
SP02025	25 mL	2/10 mL		Plastic Bags, 25 pcs/bag, 10 bags/case

50mL Serological Pipettes (Sterilization)



Biosima Model	Volume	Scale	Color Ring	Packing
SP01050	50 mL	1/2 mL		Paper-Plastic Bags, 25 pcs/bag, 8 bags/case
SP02050	50 mL	1/2 mL		Plastic Bags, 25 pcs/bag, 8 bags/case

100mL Serological Pipettes (Sterilization)



Biosima Model	Volume	Scale	Color Ring	Packing
SP01100	100 mL	1 mL		Paper-Plastic Bags, 20 pcs/bag, 6 bags/case
SP02100	100 mL	1 mL		Plastic Bags, 20 pcs/bag, 6 bags/case

Compatible with a variety of pipettors, suitable for pipetting micro reagent from deep vessels.

Cryogenic Tube



USP Raw Material

USP grade VI polystyrene (PS) raw material.



High Temperature Sterilization

Autoclavable at 121°C/15psi.



Deep Freezing

Can be frozen and refrigerated up to -196 ° C.



Low Retention

Smooth tips, no bending, low adsorption, low residue.



Chemical Resistance

Resistant to chemical reagent during use and storage.



Low Heavy Metal Content

High quality and low heavy metal content.



Traceable

All product package is coded for traceability.

The Application Of Cryogenic Tube

- > Storage of microorganisms
- > Storage of primary cells
- > Storage of cell lines
- > Storage of blood and serum samples
- > Sample transport



Cryogenic tube frozen storage tube, three in one, boxed SBS.

What does it mean to store in a cold environment ?



Storage at low temperatures is defined as temperatures below -130°C (approx. < 140K). This is a lower temperature than the temperature at which water becomes a gas. The ice no longer recrystallizes and therefore the ice crystals do not increase in size. This ensures that the possibility of chemical reactions in the sample is minimized and morphological changes (e.g., crystal formation) are avoided.

When samples are stored in the liquid nitrogen gas phase, the liquid nitrogen vapor in the sample storage creates an inert gaseous environment that prevents changes (oxidation processes) in the sample due to oxygen in the air. Cryopreservation suspends almost all chemical reactions and is a very important step when it comes to long-term storage and preventing sample degradation. The most common method is to store samples in the gas phase in tanks of liquid nitrogen or stored in a refrigerated freezer. Suitable plastic materials and accurate thread design enable these containers to have a perfect seal and reduce the risk of sample contamination.

Storage Down To -20°C	Storage Down To -80°C	Down To -196°C Storage
- Short-Term Storage - Mid-Term Storage	- Long-term storage of non-critical samples - Mid-term storage of sensitive samples	- Long-term storage - Sample Bank - Storage cells

Screw Tube





0.5mL Screw Tube (High Cover)	Biosima Model	Volume	Tube	Cap	Sterile	Packing
 Dimension Unit: mm	SCT-050-SS-H (SCO-C-H)	0.5mL	Transparent	Transparent		500 pcs/pack,10 packs/carton
	SCT-050-SS-H (SCO-CR-H)	0.5mL	Transparent	Transparent Red		500 pcs/pack,10 packs/carton
	SCT-050-SS-H (SCO-CB-H)	0.5mL	Transparent	Transparent Blue		500 pcs/pack,10 packs/carton
	SCT-050-SS-H (SCO-CY-H)	0.5mL	Transparent	Transparent Yellow		500 pcs/pack,10 packs/carton
	SCT-050-SS-H (SCO-CG-H)	0.5mL	Transparent	Transparent Green		500 pcs/pack,10 packs/carton
	SCT-050-SS-H (SCO-W-H)	0.5mL	Transparent	White		500 pcs/pack,10 packs/carton
	SCT-050-SS-H (SCO-R-H)	0.5mL	Transparent	Red		500 pcs/pack,10 packs/carton
	SCT-050-SS-H (SCO-G-H)	0.5mL	Transparent	Green		500 pcs/pack,10 packs/carton
	SCT-050-SS-H (SCO-B-H)	0.5mL	Transparent	Blue		500 pcs/pack,10 packs/carton
	SCT-050-SS-H (SCO-Y-H)	0.5mL	Transparent	Yellow		500 pcs/pack,10 packs/carton
	SCT-050-SS-H (SCO-V-H)	0.5mL	Transparent	Purple		500 pcs/pack,10 packs/carton
	SCT-050-SS-H (SCO-O-H)	0.5mL	Transparent	Orange		500 pcs/pack,10 packs/carton
	SCT-050-SS-XH (SCO-BR-H)	0.5mL	Dark Brown	Dark Brown		500 pcs/pack,10 packs/carton
0.5mL Screw Tube (High Cover)	Biosima Model	Volume	Tube	Cap	Sterile	Packing
 Dimension Unit: mm	SCT-050-SS-HS (SCO-C-HS)	0.5mL	Transparent	Transparent	✓	500 pcs/pack,10 packs/carton
	SCT-050-SS-HS (SCO-CR-HS)	0.5mL	Transparent	Transparent Red	✓	500 pcs/pack,10 packs/carton
	SCT-050-SS-HS (SCO-CB-HS)	0.5mL	Transparent	Transparent Blue	✓	500 pcs/pack,10 packs/carton
	SCT-050-SS-HS (SCO-CY-HS)	0.5mL	Transparent	Transparent Yellow	✓	500 pcs/pack,10 packs/carton
	SCT-050-SS-HS (SCO-CG-HS)	0.5mL	Transparent	Transparent Green	✓	500 pcs/pack,10 packs/carton
	SCT-050-SS-HS (SCO-W-HS)	0.5mL	Transparent	White	✓	500 pcs/pack,10 packs/carton
	SCT-050-SS-HS (SCO-R-HS)	0.5mL	Transparent	Red	✓	500 pcs/pack,10 packs/carton
	SCT-050-SS-HS (SCO-G-HS)	0.5mL	Transparent	Green	✓	500 pcs/pack,10 packs/carton
	SCT-050-SS-HS (SCO-B-HS)	0.5mL	Transparent	Blue	✓	500 pcs/pack,10 packs/carton
	SCT-050-SS-HS (SCO-Y-HS)	0.5mL	Transparent	Yellow	✓	500 pcs/pack,10 packs/carton
	SCT-050-SS-HS (SCO-V-HS)	0.5mL	Transparent	Purple	✓	500 pcs/pack,10 packs/carton
	SCT-050-SS-HS (SCO-O-HS)	0.5mL	Transparent	Orange	✓	500 pcs/pack,10 packs/carton
	SCT-050-SS-XH-S (SCO-BR-HS)	0.5mL	Dark Brown	Dark Brown	✓	500 pcs/pack,10 packs/carton

Multiple specifications and types to choose from, meeting various needs.

Transparent/Transparent Red/Transparent Blue/Transparent Yellow/Transparent Green



○ Solid White ● Solid Red ● Solid Green ● Solid Blue ● Solid Yellow ● Solid Purple ● Solid Orange ● Dark Brown

1.5mL Screw Tube (High Cover)	Biosima Model	Volume	Tube	Cap	Sterile	Packing
	SCT-150-SS-H (SCO-C-H)	1.5mL	Transparent	Transparent		500 pcs/pack,10 packs/carton
	SCT-150-SS-H (SCO-CR-H)	1.5mL	Transparent	Transparent Red		500 pcs/pack,10 packs/carton
	SCT-150-SS-H (SCO-CB-H)	1.5mL	Transparent	Transparent Blue		500 pcs/pack,10 packs/carton
	SCT-150-SS-H (SCO-CB-H)	1.5mL	Transparent	Transparent Yellow		500 pcs/pack,10 packs/carton
	SCT-150-SS-H (SCO-CG-H)	1.5mL	Transparent	Transparent Green		500 pcs/pack,10 packs/carton
	SCT-150-SS-H (SCO-W-H)	1.5mL	Transparent			500 pcs/pack,10 packs/carton
	SCT-150-SS-H (SCO-R-H)	1.5mL	Transparent			500 pcs/pack,10 packs/carton
	SCT-150-SS-H (SCO-G-H)	1.5mL	Transparent			500 pcs/pack,10 packs/carton
	SCT-150-SS-H (SCO-B-H)	1.5mL	Transparent			500 pcs/pack,10 packs/carton
	SCT-150-SS-H (SCO-Y-H)	1.5mL	Transparent			500 pcs/pack,10 packs/carton
	SCT-150-SS-H (SCO-V-H)	1.5mL	Transparent			500 pcs/pack,10 packs/carton
	SCT-150-SS-H (SCO-O-H)	1.5mL	Transparent			500 pcs/pack,10 packs/carton
	SCT-150-SS-XH (SCO-BR-H)	1.5mL	Dark Brown			500 pcs/pack,10 packs/carton
	Dimension Unit: mm					
1.5mL Screw Tube (High Cover)	Biosima Model	Volume	Tube	Cap	Sterile	Packing
	SCT-150-SS-HS (SCO-C-HS)	1.5mL	Transparent	Transparent	✓	500 pcs/pack,10 packs/carton
	SCT-150-SS-HS (SCO-CR-HS)	1.5mL	Transparent	Transparent Red	✓	500 pcs/pack,10 packs/carton
	SCT-150-SS-HS (SCO-CB-HS)	1.5mL	Transparent	Transparent Blue	✓	500 pcs/pack,10 packs/carton
	SCT-150-SS-HS (SCO-CB-HS)	1.5mL	Transparent	Transparent Yellow	✓	500 pcs/pack,10 packs/carton
	SCT-150-SS-HS (SCO-CG-HS)	1.5mL	Transparent	Transparent Green	✓	500 pcs/pack,10 packs/carton
	SCT-150-SS-HS (SCO-W-HS)	1.5mL	Transparent		✓	500 pcs/pack,10 packs/carton
	SCT-150-SS-HS (SCO-R-HS)	1.5mL	Transparent		✓	500 pcs/pack,10 packs/carton
	SCT-150-SS-HS (SCO-G-HS)	1.5mL	Transparent		✓	500 pcs/pack,10 packs/carton
	SCT-150-SS-HS (SCO-B-HS)	1.5mL	Transparent		✓	500 pcs/pack,10 packs/carton
	SCT-150-SS-HS (SCO-Y-HS)	1.5mL	Transparent		✓	500 pcs/pack,10 packs/carton
	SCT-150-SS-HS (SCO-V-HS)	1.5mL	Transparent		✓	500 pcs/pack,10 packs/carton
	SCT-150-SS-HS (SCO-O-HS)	1.5mL	Transparent		✓	500 pcs/pack,10 packs/carton
	SCT-150-SS-XH-S (SCO-BR-HS)	1.5mL	Dark Brown		✓	500 pcs/pack,10 packs/carton
	Dimension Unit: mm					

Multiple specifications and types to choose from, meeting various needs.

Transparent/Transparent Red/Transparent Blue/Transparent Yellow/Transparent Green




○ Solid White ● Solid Red ● Solid Green ● Solid Blue ● Solid Yellow ● Solid Purple ● Solid Orange ● Dark Brown

2.0mL Screw Tube (High Cover)	Biosima Model	Volume	Tube	Cap	Sterile	Packing
	SCT-200-SS-H (SCO-C-H)	2.0mL	Transparent	Transparent		500 pcs/pack,10 packs/carton
	SCT-200-SS-H (SCO-CR-H)	2.0mL	Transparent	Transparent Red		500 pcs/pack,10 packs/carton
	SCT-200-SS-H (SCO-CB-H)	2.0mL	Transparent	Transparent Blue		500 pcs/pack,10 packs/carton
	SCT-200-SS-H (SCO-CB-H)	2.0mL	Transparent	Transparent Yellow		500 pcs/pack,10 packs/carton
	SCT-200-SS-H (SCO-CG-H)	2.0mL	Transparent	Transparent Green		500 pcs/pack,10 packs/carton
	SCT-200-SS-H (SCO-W-H)	2.0mL	Transparent			500 pcs/pack,10 packs/carton
	SCT-200-SS-H (SCO-R-H)	2.0mL	Transparent			500 pcs/pack,10 packs/carton
	SCT-200-SS-H (SCO-G-H)	2.0mL	Transparent			500 pcs/pack,10 packs/carton
	SCT-200-SS-H (SCO-B-H)	2.0mL	Transparent			500 pcs/pack,10 packs/carton
	SCT-200-SS-H (SCO-Y-H)	2.0mL	Transparent			500 pcs/pack,10 packs/carton
	SCT-200-SS-H (SCO-V-H)	2.0mL	Transparent			500 pcs/pack,10 packs/carton
	SCT-200-SS-H (SCO-O-H)	2.0mL	Transparent			500 pcs/pack,10 packs/carton
	SCT-200-SS-XH (SCO-BR-H)	2.0mL	Dark Brown			500 pcs/pack,10 packs/carton
	Dimension Unit: mm					
2.0mL Screw Tube (High Cover)	Biosima Model	Volume	Tube	Cap	Sterile	Packing
	SCT-200-SS-HS (SCO-C-HS)	2.0mL	Transparent	Transparent	✓	500 pcs/pack,10 packs/carton
	SCT-200-SS-HS (SCO-CR-HS)	2.0mL	Transparent	Transparent Red	✓	500 pcs/pack,10 packs/carton
	SCT-200-SS-HS (SCO-CB-HS)	2.0mL	Transparent	Transparent Blue	✓	500 pcs/pack,10 packs/carton
	SCT-200-SS-HS (SCO-CB-HS)	2.0mL	Transparent	Transparent Yellow	✓	500 pcs/pack,10 packs/carton
	SCT-200-SS-HS (SCO-CG-HS)	2.0mL	Transparent	Transparent Green	✓	500 pcs/pack,10 packs/carton
	SCT-200-SS-HS (SCO-W-HS)	2.0mL	Transparent		✓	500 pcs/pack,10 packs/carton
	SCT-200-SS-HS (SCO-R-HS)	2.0mL	Transparent		✓	500 pcs/pack,10 packs/carton
	SCT-200-SS-HS (SCO-G-HS)	2.0mL	Transparent		✓	500 pcs/pack,10 packs/carton
	SCT-200-SS-HS (SCO-B-HS)	2.0mL	Transparent		✓	500 pcs/pack,10 packs/carton
	SCT-200-SS-HS (SCO-Y-HS)	2.0mL	Transparent		✓	500 pcs/pack,10 packs/carton
	SCT-200-SS-HS (SCO-V-HS)	2.0mL	Transparent		✓	500 pcs/pack,10 packs/carton
	SCT-200-SS-HS (SCO-O-HS)	2.0mL	Transparent		✓	500 pcs/pack,10 packs/carton
	SCT-200-SS-XH-S (SCO-BR-HS)	2.0mL	Dark Brown		✓	500 pcs/pack,10 packs/carton
	Dimension Unit: mm					

Multiple specifications and types to choose from, meeting various needs.

Transparent/Transparent Red/Transparent Blue/Transparent Yellow/Transparent Green



○ Solid White ● Solid Red ● Solid Green ● Solid Blue ● Solid Yellow ● Solid Purple ● Solid Orange ● Dark Brown

0.5mL Screw Tube (Low Cover)	Biosima Model	Volume	Tube	Cap	Sterile	Packing
 Dimension Unit: mm	SCT-050-SS (SCO-C)	0.5mL	Transparent	Transparent		500 pcs/pack,10 packs/carton
	SCT-050-SS (SCO-CR)	0.5mL	Transparent	Transparent Red		500 pcs/pack,10 packs/carton
	SCT-050-SS (SCO-CB)	0.5mL	Transparent	Transparent Blue		500 pcs/pack,10 packs/carton
	SCT-050-SS (SCO-CY)	0.5mL	Transparent	Transparent Yellow		500 pcs/pack,10 packs/carton
	SCT-050-SS (SCO-CG)	0.5mL	Transparent	Transparent Green		500 pcs/pack,10 packs/carton
	SCT-050-SS-S (SCO-C-S)	0.5mL	Transparent	Transparent	✓	500 pcs/pack,10 packs/carton
	SCT-050-SS-S (SCO-CR-S)	0.5mL	Transparent	Transparent Red	✓	500 pcs/pack,10 packs/carton
	SCT-050-SS-S (SCO-CB-S)	0.5mL	Transparent	Transparent Blue	✓	500 pcs/pack,10 packs/carton
	SCT-050-SS-S (SCO-CY-S)	0.5mL	Transparent	Transparent Yellow	✓	500 pcs/pack,10 packs/carton
	SCT-050-SS-S (SCO-CG-S)	0.5mL	Transparent	Transparent Green	✓	500 pcs/pack,10 packs/carton
1.5mL Screw Tube (Low Cover)	Biosima Model	Volume	Tube	Cap	Sterile	Packing
 Dimension Unit: mm	SCT-150-SS (SCO-C)	1.5mL	Transparent	Transparent		500 pcs/pack,10 packs/carton
	SCT-150-SS (SCO-CR)	1.5mL	Transparent	Transparent Red		500 pcs/pack,10 packs/carton
	SCT-150-SS (SCO-CB)	1.5mL	Transparent	Transparent Blue		500 pcs/pack,10 packs/carton
	SCT-150-SS (SCO-CY)	1.5mL	Transparent	Transparent Yellow		500 pcs/pack,10 packs/carton
	SCT-150-SS (SCO-CG)	1.5mL	Transparent	Transparent Green		500 pcs/pack,10 packs/carton
	SCT-150-SS-S (SCO-C-S)	1.5mL	Transparent	Transparent	✓	500 pcs/pack,10 packs/carton
	SCT-150-SS-S (SCO-CR-S)	1.5mL	Transparent	Transparent Red	✓	500 pcs/pack,10 packs/carton
	SCT-150-SS-S (SCO-CB-S)	1.5mL	Transparent	Transparent Blue	✓	500 pcs/pack,10 packs/carton
	SCT-150-SS-S (SCO-CY-S)	1.5mL	Transparent	Transparent Yellow	✓	500 pcs/pack,10 packs/carton
	SCT-150-SS-S (SCO-CG-S)	1.5mL	Transparent	Transparent Green	✓	500 pcs/pack,10 packs/carton
2.0mL Screw Tube (Low Cover)	Biosima Model	Volume	Tube	Cap	Sterile	Packing
 Dimension Unit: mm	SCT-200-SS (SCO-C)	2.0mL	Transparent	Transparent		500 pcs/pack,10 packs/carton
	SCT-200-SS (SCO-CR)	2.0mL	Transparent	Transparent Red		500 pcs/pack,10 packs/carton
	SCT-200-SS (SCO-CB)	2.0mL	Transparent	Transparent Blue		500 pcs/pack,10 packs/carton
	SCT-200-SS (SCO-CY)	2.0mL	Transparent	Transparent Yellow		500 pcs/pack,10 packs/carton
	SCT-200-SS (SCO-CG)	2.0mL	Transparent	Transparent Green		500 pcs/pack,10 packs/carton
	SCT-200-SS-S (SCO-C-S)	2.0mL	Transparent	Transparent	✓	500 pcs/pack,10 packs/carton
	SCT-200-SS-S (SCO-CR-S)	2.0mL	Transparent	Transparent Red	✓	500 pcs/pack,10 packs/carton
	SCT-200-SS-S (SCO-CB-S)	2.0mL	Transparent	Transparent Blue	✓	500 pcs/pack,10 packs/carton
	SCT-200-SS-S (SCO-CY-S)	2.0mL	Transparent	Transparent Yellow	✓	500 pcs/pack,10 packs/carton
	SCT-200-SS-S (SCO-CG-S)	2.0mL	Transparent	Transparent Green	✓	500 pcs/pack,10 packs/carton



Durable to -80 ℃~121 ℃, safe storage under liquid nitrogen gas phase conditions.

Cryogenic Tube



0.5mL Cryogenic Tube	Biosima Model	Volume	Cap	Sterile	Packing
 Dimension Unit: mm	SCT-050-SS-W-S	0.5mL	White	✓	50 pcs/pack,20 packs/carton
	SCT-050-SS-R-S	0.5mL	Red	✓	50 pcs/pack,20 packs/carton
	SCT-050-SS-G-S	0.5mL	Green	✓	50 pcs/pack,20 packs/carton
	SCT-050-SS-B-S	0.5mL	Blue	✓	50 pcs/pack,20 packs/carton
	SCT-050-SS-Y-S	0.5mL	Yellow	✓	50 pcs/pack,20 packs/carton
	SCT-050-SS-V-S	0.5mL	Purple	✓	50 pcs/pack,20 packs/carton
	SCT-050-SS-O-S	0.5mL	Orange	✓	50 pcs/pack,20 packs/carton
1.5mL Cryogenic Tube	Biosima Model	Volume	Cap	Sterile	Packing
 Dimension Unit: mm	SCT-150-SS-W-S	1.5mL	White	✓	50 pcs/pack,20 packs/carton
	SCT-150-SS-R-S	1.5mL	Red	✓	50 pcs/pack,20 packs/carton
	SCT-150-SS-G-S	1.5mL	Green	✓	50 pcs/pack,20 packs/carton
	SCT-150-SS-B-S	1.5mL	Blue	✓	50 pcs/pack,20 packs/carton
	SCT-150-SS-Y-S	1.5mL	Yellow	✓	50 pcs/pack,20 packs/carton
	SCT-150-SS-V-S	1.5mL	Purple	✓	50 pcs/pack,20 packs/carton
	SCT-150-SS-O-S	1.5mL	Orange	✓	50 pcs/pack,20 packs/carton

Durable to -196 ℃~121 ℃, safe storage under liquid nitrogen gas phase conditions.

2.0mL Cryogenic Tube	Biosima Model	Volume	Cap	Sterile	Packing
 Dimension Unit: mm	SCT-200-SS-W-S	2.0mL		✓	50 pcs/pack, 20 packs/carton
	SCT-200-SS-R-S	2.0mL		✓	50 pcs/pack, 20 packs/carton
	SCT-200-SS-G-S	2.0mL		✓	50 pcs/pack, 20 packs/carton
	SCT-200-SS-B-S	2.0mL		✓	50 pcs/pack, 20 packs/carton
	SCT-200-SS-Y-S	2.0mL		✓	50 pcs/pack, 20 packs/carton
	SCT-200-SS-V-S	2.0mL		✓	50 pcs/pack, 20 packs/carton
	SCT-200-SS-O-S	2.0mL		✓	50 pcs/pack, 20 packs/carton
5.0mL Cryogenic Tube	Biosima Model	Volume	Cap	Sterile	Packing
 Dimension Unit: mm	SCT-500-SS-W-S	5.0mL		✓	50 pcs/pack, 50 packs/box, 4 boxes/carton
	SCT-500-SS-R-S	5.0mL		✓	50 pcs/pack, 50 packs/box, 4 boxes/carton
	SCT-500-SS-G-S	5.0mL		✓	50 pcs/pack, 50 packs/box, 4 boxes/carton
	SCT-500-SS-B-S	5.0mL		✓	50 pcs/pack, 50 packs/box, 4 boxes/carton
	SCT-500-SS-Y-S	5.0mL		✓	50 pcs/pack, 50 packs/box, 4 boxes/carton
	SCT-500-SS-V-S	5.0mL		✓	50 pcs/pack, 50 packs/box, 4 boxes/carton
	SCT-500-SS-O-S	5.0mL		✓	50 pcs/pack, 50 packs/box, 4 boxes/carton



Durable to -196 ℃~121 ℃, safe storage under liquid nitrogen gas phase conditions.

Cryogenic Tube (SBS)

Cryogenic Tube (SBS)	Biosima Model	Volume	Cap	Sterile	Packing
	SC-050-B-RS	0.5mL		✓	96 pcs/box, 80 boxes/carton
	SC-100-B-RS	1.0mL		✓	96 pcs/box, 50 boxes/carton
	SC-190-B-RS	1.9mL		✓	48 pcs/box, 50 boxes/carton

Cryogenic tube frozen storage tube, three in one, boxed SBS.

Freezing Box

Paper Freezer Box	Biosima Model	Specifications	Adapter	Packing
 Tolerance Temperature: -196 ℃~121 ℃	FB-100-W2	133X133X52mm (100 Cells, White)	1.5mL Cryogenic Tube	5 pcs/box, 10 boxes/carton
	FB-81-W2	133X133X52mm (81 Cells, White)	2.0mL Cryogenic Tube	5 pcs/box, 10 boxes/carton
Freezing Box (PP)	Biosima Model	Specifications	Adapter	Packing
 Tolerance Temperature: -80 ℃~121 ℃	FB-PP100-B1	133X133X36mm (100 Cells, Blue)	0.5mL Cryogenic Tube	8 pcs/box, 10 boxes/carton
	FB-PP81-B2	133X133X52mm (81 Cells, Blue)	1.5mL/2mL Cryogenic Tube	5 pcs/box, 10 boxes/carton
	FB-PP100-B2	133X133X52mm (100 Cells, Blue)	1.5mL/2mL Cryogenic Tube	5 pcs/box, 10 boxes/carton
	FB-PP81-B3	133X133X95mm (81 Cells, Blue)	5mL Cryogenic Tube	3 pcs/box, 10 boxes/carton
Freezing Box (PC)	Biosima Model	Specifications	Adapter	Packing
 Tolerance Temperature: -196 ℃~121 ℃	FB-PC81-A2	133X133X52mm (81 Cells, Mixed Color)	1.5mL/2mL Cryogenic Tube	5 pcs/box, 4 boxes/carton
	FB-PC100-A2	133X133X52mm (100 Cells, Mixed Color)	1.5mL/2mL Cryogenic Tube	5 pcs/box, 4 boxes/carton
	FB-PC81-A3	133X133X95mm (81 Cells, Mixed Color)	5mL Cryogenic Tube	5 pcs/box, 4 boxes/carton

Multiple specifications and types to choose from, meeting various needs.     

Customized



Uniqueness
Exclusive agency agreement.

Adjustable
Adjust the parameter according to customer's requirements.

Cost-Effective
Replacement of imported consumables.

Stable Supply
Short delivery time, no stocking pressure.

DIY
Different mold and size process according to specific needs.



Customization - the key to industry upgrading.



Bize Ulařın ;

0542 246 74 62

**Herhangi bir yardım
iin bize mail atın:**

info@bio-sima.com

